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### VAROVÁNÍ

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## STOMATOLOGY ROTARY INSTRUMENTS

### CAUTION

- The caution is valid for these instruments: carbide burs, finishing burs, diamond instruments, steel and carbide cutters, instruments for root canals treatment, aluminium stands and boxes.
- Don't use agents containing phenol and hydrogen peroxide  $H_2O_2$  for disinfection.
- Instruments made of NiTi can be disinfected with 5 % solution of sodium hypochlorite NaClO, however not longer than 5 minutes.
- Attention, the endodontic instruments are sharp!
- Avoid of using of alkalic solutions (pH higher than 8), acidic (pH lower than 6), sodium hydroxide solutions or solutions containing mercury salts for instruments and boxes made of aluminium and aluminium alloys.
- The instruments made of carbide mustn't be washed in washing machines.
- Avoid of the alkalic solutions and hydroxide peroxide solutions and using of aerate water by washing of instruments with carbide parts.
- Avoid of using of high-alkalic and high-oxidative agents (peracetic acid, sodium hypochlorite) on the instruments with diamond surface.
- Remove the clinging rests of tissues with compressed air, rinse in fully desalted or distilled water.
- Used instruments for root canals treatment can cause a package perforation; this is the reason why they must be put in safety solid packages – PET bottles, canisters etc. during the waste classification.
- Only the dentists can work with the instruments.
- The producer is not responsible for injuries either on health or possessions of the patient.
- The instruments contain nickel that can cause an allergic reaction to the patient.

### Limit of repeated using:

The end of lifetime is determined through their wear and damage during using. Instruments Barbed Broaches are disposable. The manufacturer does not recommend repeated using of the instrument from the reason of possible contamination or infection transfer. The cleaning, that follows, cannot completely remove possible product contamination.

### INSTRUCTION

#### How to use:

#### Carbide burs, finishing burs, diamond instruments, steel and carbide cutters

- Choose the suitable shape of instrument for the treatment. The incorrectly chosen instrument causes an imperfect shape of preparation.
- Use the advised instrument revolutions. High revolutions affect the instrument stability and the vitality of tooth pulp in its nearness. The same is valid for insufficient cooling.
- Ensure a perfect water spray cooling (min. 50 ml/min) for all revolutions over 1500/min.
- The cooling of prepared surfaces can negatively influence the incorrectly adjusted suction. A heat defect of a tooth can happen as a result.
- Work with compressive force of 0.3–2 N according to the kind of preparation. Using of

- high compressive force can result in defect or destroying of the instrument.
- Instruments with small diameters wear out in the front part as soonest. Work with these instruments as gently as possible and with low compressive force. So you prevent the instrument from lifetime decrease.
- Use the whole surface of the function end for preparation with turbine instruments not to get to their damage.
- Check the instrument before and also after the use. Put immediately away all blunt, bent or in other way defected instruments.
- The instruments must be disinfected, washed, dried and sterilized before the first using and immediately after their using.

#### Instruments for root canals treatment: Barbed Broach, K-Reamer, K-File, W-File, H-File, M-Reamer, M-File, Spreader, Plugger, Smooth Broach, Paste Carrier, Rasp, Gutta Condensor

- Work carefully and follow all the rules of endodontic practice.
- Instrument M-File is only intended for smoothing with a file without rotation and that counts even in case that the instrument is grasped in the dental tenon. The instrument is unique in the market for its cut and it is also patented. Great sharpness ensures superior treatment of the root canal.
- K-File can be recommended for smoothing with a file using partial rotation.
- Instrument M-Reamer is unique in the market for its cut and it is also patented. It is necessary to work carefully with these instruments at low revolutions and to follow all the rules used in endodontic practice.
- The Gutta Condensor is used to mechanically condense gutta percha in the root canal.
- Work at low revolutions.
- Ensure that the access to the root canal would be large and smooth enough.
- Measure the working length of the canal.
- Work from small diameters to the large ones.
- Do not omit any size of the root instrument.
- The instruments must be disinfected, washed, dried and sterilized before first use and immediately after the use.

#### Instructions for disinfection, cleaning and sterilization of dental rotary instruments

Where to use them: The instructions are valid for all above mentioned instruments.

From hygienic and safety reasons must be all instruments disinfected, washed and made sterile before intra-oral use.

#### Preparation for cleaning: Any specific requirements.

**Automated cleaning:** Carefully respect the instructions and information about advised cleaning agent concentrations. Read attentively all warnings in the introduction of Instruction for use. Carefully check the instruments after finished washing cycle, like after hand washing.

**Cleaning by hand:** Clean away the rests of dentine with soft brush, rinse and dry the instruments. Clean away the possible clinging dirt with compressed air after drying.

**Disinfection:** Immerse the instruments into disinfection agent immediately after their using or as soon as possible. The disinfection agents must be used in conformity with instructions mentioned on the label. For the preparation of disinfection and washing

solutions respect the procedure given by their producer. High concentration can result in appearance of chinks or lead to a damage of coloured distinguishing rings. Respect the disinfection given time.

*Recommended agents:* Sekusept pulver, Neodisher septo MED.

**Drying:** If the drying is a part of the washing and disinfection device cycle, the temperature of 120 °C mustn't be exceeded. Protect the instruments of carbon steel with an anticorrosive agent.

**Check and function testing:** Carefully check and put away all defected, bent, broken, blunt, corroding or other destroyed instruments. Put the checked instruments into dishes or stands. Treat the instruments of carbonaceous steels with an anticorrosive agent.

**Packing:** The instruments will be put separately into packages or into sterilizing dishes that are wrapped after it. The stands with instruments sterilize separately, open the boxes and sterilize them also individually. The standard packing material can be used. Ensure that the package is big enough and the instrument can be put in it without seal stress. Write over the packages.

**Sterilization:** If it isn't mentioned in "Special warning" otherwise, we recommend the sterilization with wet heat in steam sterilizer equipped with antibacterial filter at the temperature of 121 °C, overpressure of 205 kPa and for 20 minutes or at the temperature of 134 °C, overpressure of 304 kPa for 10 minutes.

If there is advice for hot-air sterilization in "Special instructions", it is performed in devices with forced air circulation at the temperature of 160 °C for 60 minutes or at the temperature of 170 °C for 30 minutes or at the temperature of 180 °C for 20 minutes, if it isn't mentioned otherwise. The hot-air sterilizer is opened when the cycle is finished and cooled down to at least 80 °C.

Packed sterilized instruments can be sterilized repeatedly.

We recommend to sterilize the steel instruments separately from other material instruments and to protect them with anticorrosive agent.

When sterilising more instruments at one cycle, it is ensured so that the maximum load of the steriliser is not exceeded.

**Storage:** Store in dry spaces and prevent from the repeated contamination.

**Other information:** The joint-stock company MEDIN validated the above mentioned methods of cleaning, disinfection, drying and sterilization. Through the validation were all the procedures recognized as qualified for the preparation on medical instruments for repeated using.

The user is responsible for ensuring of necessary materials, equipment and staff training in accordance with these specified procedures so that the required effect could be reached.

The user can use some other methods, but the effectiveness of these alternative methods must be verified by specialized laboratory.

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