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## MILLING CUTTERS, DRILLS AND OTHER ROTARY INSTRUMENTS

## Description of instruments

## Drills

These are rotary surgical instruments used for drilling of openings into the bone and hard tissues during surgeries. The clamping part is adapted for clamping of various types of drilling-machines. Some drills are equipped with marks for orientation of the surgeon on the depth of drilling.

## Milling cutters

These are rotary surgical instruments used for milling of bones and hard tissues. The clamping part is adapted for various types of drilling-machines.

## Special milling cutters

ORL: Used for drilling and widening of an opening in the bone or other hard tissues during surgeries of the middle ear cavity or paranasal sinuses. The functional part of milling cutters has a spherical, pear like or conical shape. The clamping part is adapted for various types of drilling-machines.

Marrow milling cutter: Used for pre-drilling of a bone opening prior to the insertion of the pin. The milling cutters are flexible.

## Screw-taps

Used for making of a thread in a pre-drilled opening in the bone. The shape of the thread corresponds with that of the thread of the inserted screw (cortical, spongy etc.). The clamping part is adapted for clamping of various types of drilling-machines.

## Compatibility

The products are compatible with a drilling-machine and tangs.

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- making an opening in the bone to introduce a bolting screw
- making an opening to introduce a cannula screw

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- making a driving opening to introduce an implant or a bone graft
- making an opening into the bone during spinal surgeries (Casper)

## Milling cutters and special milling cutters:

Used in surgery and traumatology for

- making a driving opening to insert an implant or a bone graft
- acquiring a bone cylinder from the bone to implant it into a required site

## Screw taps:

Used in surgery and traumatology for

- making an opening and a thread to insert a screw

## Contraindications

No contraindications have been determined.

## Procedure of using

Fasten the instrument into an electric drilling-machine and use for a required operation, applying adequate pressure.

While using the instruments, use protective cases.

## WARNING

- Prior to its use, always check the tugging of the instrument visually.
- After the instrument becomes blunt, ensure its sharpening.
- Attention! The instrument is sharp.
- The instrument is equipped with edge protection prior to its delivery and for further handling.
- While the instrument is sterilised, take care that its edge is not damaged when in contact with other instruments; you will therefore increase the lifetime of the instrument.
- Prior to its use, always check the quality of the edge visually, not by touch.
- The instrument is indirectly connected with a source of electric energy. As the drilling-machine is on, avoid its other touching with the instrument than that for a given purpose.
- Prior to the use of the instrument, make sure that its clamping part corresponds with an opening in the drilling-machine or with a drilling-machine shaft.

## Limit of repeated treatment

If the manufacturer's recommendations on storing are followed, the lifetime is limited to 20 years due to the lifetime of containers and filing of documentation. As the instrument has been used, its lifetime is given by its wearing. Some instruments can be repaired by sharpening their edges.

## INSTRUCTIONS

**PRECAUTIONS:** While using steel instruments, it is necessary to avoid agents containing chloride ions which may cause corrosion;

**Place of use:** Professional workplaces of healthcare providers, particularly surgical and orthopaedic

**Storing and transport:** Protect the surgical instruments against deformation and abrasion; avoid their contact with other metal materials

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and chemicals while handling.

## Preparation for cleaning

During cleaning and disinfection, it is necessary to avoid agents containing chloride ions which can cause corrosion of steel parts.

1. Use mechanical cleaning with water by means of brush.
2. Sluice the instruments with pressure water.
3. Immerse the instruments into a disinfection solution for 20 to 30 minutes. The recommended disinfection agent is Sekusept pulver, NeodisherseptoMED. In the case of using another solution, there is a risk of damaging the instruments.
4. Repeatedly sluice with a stream of pressure water unless uncontaminated water is running. It is possible to make sluicing in conventionally supplied pressure washing machines after linking of the instruments to jets.

**Disinfection:** It is recommended to use either a washing device and thermic or thermochemical disinfection or hand washing of the instruments after their chemical disinfection with antiviral effect; it is possible to join both phases here, using disinfection agents with combined washing effect. The disinfection agent can be used in accordance with the instructions stated on the label. While disinfection and washing solutions are being prepared, it is necessary to comply with the procedure stated by the manufacturer. The instruments intended for dismounting are disassembled and each part is further considered an individual instrument. Recommended agents: Sekusept pulver, NeodisherseptoMED. **Drying:** If drying is a part of washing and disinfection device cycle, the temperature 120 °C must not be exceeded.

**Maintenance:** After each operation has been accomplished, it is necessary to check the instruments. In case of their damage, it is necessary to deliver the components to the manufacturer to make a repair.

**Checking and test of function:** Prior to its use, always make a visual check of tugging of the instrument.

Prior to its use, always check the quality of the edge visually, not by touch.

Prior to the use of the instrument, make sure that its clamping part corresponds with an opening in the drilling-machine or with a drilling-machine shaft.

**Sterilisation:** Unless otherwise stated, we recommend wet heat sterilisation in a steam steriliser equipped with an antibacterial filter at the temperature 121 °C and overpressure 205 kPa for 20 minutes or at the temperature 134 °C and overpressure 304 kPa for 10 minutes.

**Other information:** The joint-stock company MEDIN validated the above-mentioned procedures of cleaning, disinfection, drying and sterilisation. By the validation, these procedures were recognised as eligible for the preparation of medical devices for repeated use.

The user is responsible for securing of necessary materials, equipment and staff training in accordance with these specified procedures to reach the required effect.

The user can also use other methods but the efficacy of these alternative procedures must be verified by a professional laboratory.

**Storing:** Packed products are stored in dry and clean environment without extreme temperatures and humidity, out of reach of direct sunlight.

**Disposal:** While being disposed, the used instruments can cause perforation of the packing; while being sorted out, they must therefore be put in safe solid containers – PET bottles, small tanks etc.

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