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## LIPOSUCTION CANNULAS, LIPOSUCTION AND FAT TRANSPLANTATION WITH SYRINGE, ACCESSORY

### Liposuction cannulas with fixing into the hose

The liposuction cannulas are meant for suction of fat from body through small nygma in concrete point after fixing on the hose of suction apparatus. Cannulas are pipes of various length and diameter with blunt full spike on the end and with one or more gaps one after another (or around). On the opposite end is the cannula equipped with handle with knurled part for fixing on the hose. There is the surface for polex holding the cannula – the surface position is opposite to the cannula gap. With moving of cannula in the body backwards and forwards is the fat abraded and with suction taken into suction hose.

### Liposuction cannulas with fixing into syringe

The cannulas in their function part correspond in the shape to cannulas with handle (of smaller sizes) but the fixing is adjusted to syringe of Luer LOCK type with volume of 50–60 ml. For using of suction with syringe secure the attached cannula with handle on the syringe. The suction in syringe secures the forceps which is attached on the opposite side of Luer LOCK type syringe with volume of 50–60 ml. By twist of syringe piston the sharp edge of forceps cuts into piston and so secure its position. After making of low pressure in the syringe will it be able to work with cannula.

### Liposuction fat transplantation with syringe

For implantation of own patient's fat will be sucked small fat volume with blunt cannula (the best of diameter 3 mm, with one gap in length of 150 mm) from small nygma in adequate size of cannula diameter. This material will be bypassed into some other syringe and several times flushed through bypassing with physiological solution to put aside the liquids and to retain only stroma with fat stored in cells. After the last douche the syringe will be put on the stand and the taken fat will be let to settle. After fat sedimentation and squirting off of free liquid the necessary fat volume will be applied with implant needle.

The forceps secures the low pressure in syringe. The holder (teflon collar) of syringe spike and cannula base in the connection place are used for fixing of the connection.

Cleaning of cannulas after their using is done with mandrels of diameter adequate to the cannula gap.

### Accessory

- 1) the cannula handle – teflon collar fixing the coupling of cannula and syringe
- 2) securing forceps – on the syringe of Luer LOCK type with 50–60 ml of volume – to order at the manufacturer of cannulas (description of using see above)
- 3) bypassing connecting piece – small tube with ends of connection of Luer LOCK 50–60 ml type, for sterile bypassing of aspirated material into some other syringe
- 4) stand – with fixing to the syringe of Luer LOCK 50–60 ml type is used for putting away of full syringe into vertical position for fat sedimentation

Notice: The choice of cannulas is available in special catalogue at the manufacturer.

## LIPOSUKČNÉ KANYLY, LIPOSUKCIA A TRANSPLANTÁCIA TUKU INJEKČNOU STRIEKAČKOU, PŘÍSLUŠENSTVO

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