VAROVÁNÍ

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CAUTION

PRECAUTION:

- The implant must be never used again! The previous using can decrease its fatigue strength and when used again the osteosynthesis can fail. The used implant can be contaminated and its repeated using can cause inflammatory reaction of the organism. It can also transfer the illness that the first patient suffered from.
- The implant cannot be unlimitedly resistant to the level of functions which are in accordance with normal healthy bone. It can be used as long as the bone is not healed, the implant for fracture fixation at active younger patients is usually removed after this healing time. However at the older patient the surgeon can consider the implant leaving for fracture fixation after the healing than risk the additional surgical intervention of its evulsion even if the implant loses its function yet.
- The surgeon must be acquainted with the operation process by using of specified instruments.
- Use only the specified instruments for surgical procedure; it mustn't be injured to be used well.
- Consider the use of implants for patients where you aren't sure that they will cooperate with doctor after operation (e.g. drunkards, drug-addicts, mental sick).
- Protect the implant from deformation and scratching, avoid the contact with other metal materials and chemicals by manipulation.
- If the adjustment of implant shape is necessary, make it only with specified way.
- Don't use various material-different implants for one patient, because the contact of two various metals could speed up the process of corrosion what could result in earlier drop of the implant effect and lifetime, further in increased volume of metal alloys that release in the patient's body and could cause a toxic reaction.
- Decide the possible untimely implant removing if you surmise the bad bone healing or that the implant can lose, bend, break in the bone, because so can it get to heavy patient's injury.
- Dispose of all shapes and sizes of implants in the medical workplace which use can be supposed.
- Take into consideration the temporary preventive means that reduce the weight loading for the purpose of minimisation of early loading of the mean for fracture fixations as long as the evident bone junction isn't possible to see.
- Decide the implant removing for internal fixation with taking account of factors that are risks for patients in consequence of additional surgical procedure and troubles during and after implant removing.
- Patients which smoke or usually use product with nicotine have escalated number of uncurd fractures. Patients like this must be warned of this fact and eventual consequences. Occurrence of implant failure is rather probable at heavy, physically active person.
- The producer does not bear the responsibility if the product is used in a different way than it is mentioned in the Instructions for use or in the Surgery procedure.
- The producer does not bear the responsibility, if the product is used in combination with implants of other producers.

Information for patient:

- Warn the patient that he must exactly follow the instructions of nursing doctor.
- Instruct the patient of operation process, possible risk and rehabilitation program.
- Instruct the patient when he can fully load the bone with internal fixation in order to avoid the implant damage.
- Warn the patient that the bone is predisposed to injury during treatment above all in points when the implant ends.
- Announce to the patient when the implant is supposed to be removed.
- Instruct the patient that he should abide adequate activity level after implant removing because it can get to repeated fracture before the bone is wholly healed and the bone cavity filled.

Possible side effects:

- Non-executed junction of fracture ends, osteoporosis, repressed revascularisation and bad ectostosis result in loosening, bending and breaking of implant or early loss of firm bone fixation.
- Bad junction of fracture point following from wrong equalization.
- Heightened reaction of fibrous tissue around the fracture point because of unstable comminuted fractures.
- Early or late infection, both deep and surface.
- Avascular necrosteon.
- Shortening of affected bone side (fracture)
- Subclinical injury of nerves can appear as a result of surgical trauma.
- Allergy to Ni.
- Risk of faulty diagnosis by dint of diagnostic devices CT or MRI where the occurrence of steel implant in the given locality causes wrong view.

The heavy complications can occur with every surgical operation for aftermath elimination of multiple injuries. These complications include except of that also urogenital, gastrointestinal and vascular injuries including thrombosis, bronchopulmonary injuries including embolism, heart failure and death.

Limit of repeated treatment: The implants are meant for humane disposable use. In case, the instruments were ready only for a surgical operation and there were not in direct contact with the patient, they can be prepared again for using by the process of cleaning, disinfection and sterilization – described further. Implants can be repeatedly disinfected and re-sterilized.

INDICATION

The implant is meant for holding of fractured bone surface in mutual contact for the purpose of fracture healing facilitation. The indications of concrete implants are mentioned on relevant separate package leaflet.

CONTRAINDICATION

1. Deficient quantity or quality of bone that can hinder the correct bone fixation.
2. Whichever fully developed or supposed latent infection.
3. Previous course of infections.
4. Whichever patient's mental disorder or neuromuscular disease that would lead to unacceptable risk in connection with fixation failure or with complications in the postoperative phase
5. Aggravated vascularity that would hinder the necessary blood distribution to the fracture or operating point.
6. Malignancy in the fracture point.
7. In all these cases when the instrument could come into conflict with anatomical or physiological structures.
8. Use of steel implant if the patient is allergic to Ni.

INSTRUCTIONS

Atention: The implant is packed in polyethylene closed packaging. It is cleaned and disinfected from the producer. Take the implant out of the packaging before first using and sterilize it according to the directive. If the packaging is damaged, it is necessary to take the implant back to the producer for repacking.

Where to use it: Professional workplaces of health care providers, mainly surgical and orthopaedic.

Storage and transport: Protect the implant from deformation and scratching; avoid the contact with other metal materials and chemicals during manipulation.

Preparation for cleaning: The implants are cleaned by the producer and therefore there is no need to clean them or disinfect. If the disinfection is still needed on unused implant, it is necessary to use such disinfectant that does not contain any chloride ions and is suitable for steel implants. The implant must not be damaged (scratched) during disinfection.

Disinfection: It is recommended either the device for washing and for thermic or thermochemical disinfection or washing of implants by hand after their chemical disinfection with some agent of virucide effect, in the process can be both phases made together by using of disinfection agents with combined washing effect. The disinfection agent can be used in conformity of instructions mentioned on the plate. The procedure for preparation of disinfection and washing solutions given by producer must be respected. Implants meant for disassembly will be de dismount and each part will be further taken for separate implant. When disinfection is finished, the implants are sluiced under non-mineral water, for eliminate of chemical residues from all surfaces.

Recommended agents: Sekusept pulver, NeodisherseptMED.

Drying: If the drying is a part of washing and disinfection device cycle, the temperature of 120 °C mustn't be exceeded.

Maintenance: Without maintenance.

Control and function check up: If the cleaning and disinfection is performed before usage of the implant, it is necessary to check visually whether the implant was not damaged.

Packing: The standard packing material can be used. Secure that the package is big enough and the implant can be put in it without seal stress. The implants will be put separately into packages or into sterilizing dishes that are wrapped after it (cutting edges must be protected).

Sterilization: We recommend to sterilize the implants separately of other instruments.

If it isn't mentioned in Special instruction otherwise, we recommend the sterilization with wet heat in steam sterilizer equipped with antibacterial filter at the temperature of 121 °C, over-pressure of 205 kPa for 20 minutes or at the temperature of 134 °C, overpressure of 304 kPa for 10 minutes. If an advice for hot-air sterilization is mentioned in Special instructions, it is provided in devices with forced air circulation at the temperature of 160 °C for 60 minutes or at the temperature of 170 °C for 30 minutes or at the temperature of 180 °C for 20 minutes, if it isn't mentioned otherwise. The hot-air sterilizer will be opened later on cooling down to 80 °C after the cycle finishing.

Storage: Store the implant in dry spaces without exposing to direct effect of chemicals.

Next information: For the sterilization of more implants in one autoclave cycle must it be guaranteed that the maximal volume of sterilizer won't be overfilled. To store, transfer and sterilize medical means it is possible to use cases, boxes, sieves and stands. Imperforated cases, boxes and sieves are not possible to be used as containers for sterilization of medical means. Open them before sterilization.

Liquidation of implant: After use is implant put down as a dangerous waste. The user is accountable for implementation of procedure of harmless manipulation and liquidation of product. After dry the used implants are remove as a potential dangerous waste.

The joint-stock company MEDIN validated the above mentioned methods of cleaning, disinfection, drying and sterilization. Through the validation were all procedures recognized as QUALIFIED for the preparation on medical devices for repeated using. The user is responsible for securing of necessary materials, equipment and staff training in accordance with these specified procedures to reach the required effect. The user can use some other methods but he must consider and judge (evaluate) the effectivity of these procedures considering the potential unwanted impacts.

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