CZ

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KIRSCHNEROVY DRÁTY, K-DRÁTY

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EN

Information on use PL0109
READ CAREFULLY!

KIRSCHNER WIRES, K-WIRES

Kirschner and K-wires are stiff metallic wires that are used to achieve the inner traction or fixation of fracture, to stabilize the bone fragments; they are also used as guide wires for pre-reaming holes and cutting threads in bones, for following use of cannulated bone screws and dynamic compression screws. The wires with threads are possible to use to capture and attract the bone fragments.

INDICATIONS:

- inner traction or fixation of bone fractures, stabilization of bone fragments
- osteotomy
- arthrodesis
- deformity correction
- revision procedures, when other method were not successful
- guide wires for pre-reaming holes and cutting threads
- guide wires of cannulated bone screws and dynamic compression screws (smooth wires)

CONTRAINDICATIONS:

The operating surgeon must be well acquainted with the technique of operation.

Conditions that increase the risk of failure:

- active or preconditioned latent infections or inflammation in operative area or close to it
- sensitivity to implant material (documented or assumed)
- obesity or overweight of the patient
- mental or neuromuscular disorders of the patient that may create unpredictable risk of fixation, especially in postoperative period
- diminished vasculature that may decrease proper blood supply of the area operation

Recommended procedure operation

We recommend introducing the wires using a reamer. The flat fastening wire ends help better transmission of twisting moment. K-wires have better reaming capacity than Kirschner wires and they can be used for pre-reaming of holes instead of bone reamers.

Before using, it is necessary to check whether the wire is not bent or damaged otherwise. During introduction by a drill, the wire must be pushed in continuously, equally and in a linear way.

Manufacturer: CZ43378030 *Number and date of last revision: R02/2011-11-30*
MEDIN, a.s. | Vlachovická 619 | CZ 592 31 Nové Město na Moravě | Czech Republic  0434
tel.: +420 566 684 336 | fax: +420 566 684 385 | e-mail: export@medin.cz

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