CZ

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**EN**

Information on use PL0112
READ CAREFULLY!

SCREWDRIVER SAFETY AND PULL

PURPOSE

These are two instruments, which complement the instruments for nailing of the lower limb. They serve for implementation or extraction of safety screws 129 79 4670÷4800 and 129 79 4810÷4940.

Indications and contraindications:

They depend on the respective type of locking screws.

INSTRUCTIONS FOR USE

Product Use

Safety screwdriver (further only SSD) is one instrument, which enables not only deduction of the twisting moment but also the deduction of the axial force at the implementation of the safety screws. SSD has to be well reassembled before using. Out of standard implementation the pipe has to be hold for spot with the opened key of the size 10 and the matrix has to be tightened by the hand in the clockwise.

The clamp of the safety screw is provided so, that the hexagon is put into the head of the safety screw and with the hand rotation of the matrix (the pipe is also rotated) is the matrix easy tightened.

ATTENTION: There is the left-hand thread on the end of the pipe, what mean, that the safety screw is clamped through the rotation of the matrix counter clockwise and loosed with the clockwise rotation. Please make the clamping and loosening only with the hand not with the key.

Operator implements so clamped safety screw on the SSD into the competent pre-drilled hole without deducing of other force than the axial force and the twisting moment. The hole pre-drilling for safety screws 5 is provided with the drill of the 5 mm diameter through both corticalis. After this will be implemented the loosening of the safety screw with the clockwise rotation of the matrix.

Before washing, disinfection and sterilisation is good to make the safety screw and so, that the pipe will be hold with the key of the size 10 and by the hand counter clockwise rotation the pipe will be loosed.

Pull screwdriver is an instrument, which enable the extraction of above mentioned safety screws. The using is very simply. The head of the extracted screw is the best cleaned away the soft parts. The pull screwdriver will be put on the axis of the safety crew. With its counter clockwise rotation the end of the screwdriver will be turned into the head of the crew and with next rotation the safety screw starts again to loose.

It is possible to work on the screwdriver with the axial force and so more easily delete the safety screw.

INSTRUCTIONS

Place of use: Specialized facility for health-care providers.

The procedures of cleaning, disinfection, drying and sterilization can be found in general instructions PL0169.

Manufacturer: CZ43378030 Number and date of last revision: R01/2012-01-24
MEDIN, a.s. | Vlachovická 619 | CZ 592 31 Nové Město na Moravě | Czech Republic
tel.: +420 566 684 336 | fax: +420 566 684 385 | e-mail: export@medin.cz

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