CZ

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ÚČEL POUŽITÍ

Popis:

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PURPOSE

Description:

Retentive cups and reticles are used in alloplasty of hip joint when it is not possible to anchor a plastic plate using cement when bone is missing or it is damaged in area of acetabulum.

Indications:

Used in cases of:
– replacement of pelvis
– osteoporosis
– acetabular protrusion
– joint stiffness
– old pelvis structure

Conditions for choice of a patient:

– relieve pain and improve function
– adaptability of a patient and his/her goodwill to follow the medical instructions, among others weight control and physical activity
– balanced diet of the patient
– patient’s skeleton must fully developed

Contraindications:

Absolute contraindications: infection, sepsis, osteomyelitis.

Relative contraindications:

– non-cooperating patient or patient with neurological problems, e.g. unable to follow instructions
– severe osteoporosis
– metabolic problems disturbing the bone structure
– osteomalacia
– distant inflammation focuses that may invade the area of implant
– serious joint destruction, distinct bone loss or resorption of bone that is visible on x-ray image
– vascular insufficiency
– muscular atrophy
– neuromuscular diseases

Conditions that increase the risk of failure:

– active or preconditioned latent infections or inflammation in operation area or close to it
– diminished vasculature that may decrease proper blood supply of area of operation
– diminished bone tissue disturbed by a disease, infection or previous implantation that cannot support the implants
– sensitivity to implant material, documented or assumed
– obesity or overweight of the patient
– non-standard anatomy of the patient in the operation area, especially bone anatomy
– mental or neuromuscular disorders of the patient that may create unpredictable risk of fixation, especially in postoperative period

Operating technique:

The operating surgeon must be well acquainted with the technique of operation.

Surface unevenness shall be treated by cutter. Cup or reticle is wedged firmly in acetabulum so that the hemisphere is in close contact with medial wall of acetabulum. It can be achieved by choosing the right dimensions. In protrusion it is usually necessary to perform following osteoplasty. The cups and reticles are fixed by screws. In hard bone cortical bone screws shall be used, either with thread precutting of 5 to 10 mm or without precutting. In soft bone cancellous bone screws shall be used directly (the best implementation with help of screwdriver with cardan joint). When the reticle is used, the first hole is made in its mid point. Plastic joint shell which is the same size as cup or reticle is first settled for a test. That is followed by actual implantation of joint shell.

It is always essential to use special instruments for the operation.

Notes:

– recommended bone cements: CMW manufacturer De Puy, APP manufacturer Alinex bone cements and pharmaceuticals company Biomet. Company MEDIN validated these medical devices for their products.

WARNINGS:

When using other than bone cement company MEDIN recommended before surgery is necessary to ensure compatibility and retention of bone cement basket or nets which are made of implant steel according to ISO 5832-1, or including the acetabular cup. Then we recommend to verify the compatibility in the package insert bone cement or check with the supplier of these components.

Incorrect tailoring, positioning, orientation and fixation of implants or non-standard load may decrease fatigue life of implants. Incorrect preoperative and preoperative manipulation with the implant and its damage (scratches and notches) may lead to corrosion and fatigue failure. It is essential to use clean gloves to manipulate with the implants. The operating surgeon must be well experienced in the technique of operation. It is forbidden to model the implants. Bending may reduce the strength of metal. If there is material discontinuity in the area of implantation insufficient support of the implants may lead to failure of implants or the whole procedure. Screws must be fully shouldered to secure a stable fixation. Complete perforation of pelvic bone must be excluded in order not to endanger the bodily structures inside pelvis.

Possible adverse effects: sensitivity to implant material, early or late postoperative infection or allergy, bone perforation, implant migration or loss, periarticular calcification or ossification, undesirable limb shortening, intraoperative or postoperative fractures and pain.

INSTRUCTIONS

The procedures of cleaning, disinfection, drying and sterilization can be found in general instructions PLO088.

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