CZ

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tel.: 566 684 327 | fax: 566 684 384 | e-mail: prodej@medin.cz | www.medin.cz

EN

Information on use PL0139 READ CAREFULLY!

BONE SCREWS, WASHERS

Identification HA – cortical screws

Identification HB – cancellous screws

Choice of sizes, relevant opening sizes, instruments for opening according to the catalogue of MEDIN, a.s.

Bone screws (type ball with inner hexagon or a cross groove) HA 2; HA 2.7; HA 3.5; HB 4; HA 4.5; HB 6.5

The screws are used in traumatology and orthopaedics for immobilisation of bone plates at the bone, for fixing of nails.

Material: implant steel or titanium

a/ not-self-tapping version

b/ self-tapping version

Canalized bone screws (type ball with inner hexagon) HB 4.5; HA 5; HB 7

The axis of the screws is made with a longitudinal opening for the line of extension wires; along the wires inserted into the bone during the operation are lead the instruments for windowing and threading for following screwing in.

Indication: the fractures of calf bone, internal ankle joint, calcaneus, distal humerus, ulna, radius, knee-pan, distal tibia, foot root bone, fixation of ligaments, circular pelvis fracture, intrathecal fracture of femoral jugulum, intratrochanteric fracture of upper meros, tibia jugulum fractures, divulsion of ileosacral joint, astragalus joint arthrodesis.

Material: implant steel.

Maleolar screws (type ball with inner hexagon, thread on the screw-half) HA 4.5

The screws are used in traumatology and orthopaedics for fractures of internal ankle joint, calcaneus, tibia jugulum and the like.

Canalized bone screws – special use (without head, with inner hexagon) HB 6; HB 7; HB 8; HB 9

The screws are used for fixation of graft from patella during the reconstruction of decussate knee ligaments.

Material: titanium. The screws are part of the set – see the catalogue of MEDIN, a.s. – using instruction for the set no. PL0115 on request at the manufacturer.

Canalized bone screws – special use (without head, with inner hexagon) HB 7

The screws are used for fixation of graft during the reconstruction of ACL with using of semitendionus ligament graft; for the screw fixation in tibia will be used a special fixating washer – clip.

Material: titanium. The screws are part of the set – see the catalogue of MEDIN, a.s. – using instruction for the set no. PL0115 on request at the manufacturer.

Clips are used for laying under the screw for fixation of graft during the reconstruction of ALC

Material: titanium. The clips are part of the set – see the catalogue of MEDIN, a.s. – using instruction for the set no. PL0115 on request at the manufacturer.

Washers are used by osteosynthesis under the screw head.

Material: implant steel.

Washers under the binding wire: are used by osteosynthesis by dint of binding wire that is threaded through the washer opening and fixing to the bone together with the washer.

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