DLAHY ACETABULÁRNÍ

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ACETABULAR PLATES

Acetabular plates serve for stabilization of damaged pelvic bones. They consist of a flat strip which can be divided into three parts: an insertion one, with an opening for the pilot of the control lever; a fixation one with threaded openings for cortical screws; and a formable one with a system of lightening holes, notches and slots so that it can be easily formed according to the individual curvature of the bone. This part is pre-bent along the oblique notch, approximately at the last third of the plate's length, corresponding to the area of linea arcuata of the pelvic bone. The set of instruments for acetabular plates serves for pre-shaping and insertion of plates and for preparation of the pelvic bone for being connected to the plate by means of screws but it contains only instruments that cannot be found among common equipment of operating theatres. The aiming arm consists of an insertion handle with a groove for a plate, a control lever with a pivot for fixing of the insertion part of the plate, a sliding arm with drill sleeves for assembling of the drilling set of instruments, a screw for arresting of the sliding arm and a screw for connecting of the control lever and insertion handle. The sleeves, depth gauge, drill and screwdriver are of common shapes but of non-standard dimensions, which are adapted to the aiming arm. Acetabular tapes contain a flat strip of the same width as the plates, but longer. They have notches for easy shaping like the plates do.

Instructions for use:

The plates are inserted by means of the acetabular insertion handle, along the inner side of the pelvic bone. Prior to the insertion of the plate, it is reasonable to use acetabular tapes of the appropriate type. The acetabular tape (left short, right short, left long, right long) is to be inserted and shaped manually according to the actual curvature of the pelvic bone. This curvature is then applied to the plate. Tapes can be used in several operations although their life is shorter than that of other parts of the set of instruments. Creating the desired curvature of the plate itself (both flexion and torsion are possible) needs to be performed gradually and carefully enough so that cyclic bending can be avoided (i.e. to avoid the need for re-straightening the plate) as it reduces the life of the plate and increases the risk of damage to the plate when placing pressure on the pelvis. The aiming instrument shall be used in an unassembled state. First, the insertion handle is assembled with the plate and the control lever by means of a connecting screw, which is to be tightened slightly with a common ratchet wrench No. 13. This system is then inserted into the pelvis so that the plate is positioned along the inner side of the pelvic bone into the same place where the acetabular tape has been shaped. The arm is mounted on the insertion handle and secured by tightening the arresting screw, again with the previously used wrench No. 13 and, into the arm's sleeves, we gradually insert the outer sleeve and mandren for marking center holes for the drill, the outer sleeve and depth gauge for resetting it and both the inner and outer sleeves for guiding the drill into the marked holes in the pelvic bone. The drill is then chucked into the drill and holes are made in the pelvic bones.

Instruction of depth drilling calculation:

The depth of the drilling is measured with the help of scales. The height of the aimer's arm is set up according to a need (the best on zero position for easier counting). After insertion of the outer sleeve into the arm, we read the value on the sleeve, we deduct the value of the height of the aimer's arm from the value of the sleeve; resultant value is the drilling depth. The scale is upside down on the drill! Cortical screws with a 4.5 thread of appropriate length are driven into the openings prepared so that they run through the full length of the thread in the plate. Only then is the plate loosened from the guiding part of the aiming arm and the aiming arm is disassembled by loosening the arresting screw and unscrewing the connection screw.

Note:

It is necessary to make a preliminary estimation of sizes of implants to be used prior to the operation. Should the patient's pelvis have atypical dimensions, the required implants can be acquired by sending an order from the MEDIN company.

We recommend to put the aimer together before the operation and to go through the instruction of depth drilling calculation to prevent from accidental plate drilling during operation and possible patient's health harm.

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