CZ

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Indication of intramedullary nail (tibial or femoral)
Nail enables the synthesis of broken parts of body – tibia or femur. It is not suitable for stabilisation of proximal or distal fractures of these bones.

Indication of Femoral intramedullary retrograde nail (further retrograde nail)
The nail enables closed locked osteosynthesis; the advantage is to use it on fractures of distal femur type A and C according to AO classification. The most crucial factor for indication of osteosynthesis by antegrade or retrograde intramedullary nail is the height of the fracture line above the ceiling of intercondylar femur foss. Indication of osteosynthesis by locking antegrade nail happens when fracture line of distal fragment interferes min. 5–6 cm above ceiling of intercondylar femur foss. Osteosynthesis by locking retrograde nail happens when the fracture line of distal fragment interferes 3–4 cm above ceiling of intercondylar foss.

The list of implants and instrumentation (see the catalogue):
No 1 Tibial intramedullary nails – steel
No 2 Femoral intramedullary nails – steel
No 3 Tibial intramedullary nails – titanium
No 4 Femoral intramedullary nails – titanium
No 5 Femoral intramedullary retrograde nails
No 6 Instrumentation

1) MEDIN Femoral intramedullary nails are hollow without groove on the surface, have a circular section. MEDIN has at its disposal nails of diameter 10–15 mm, sizes 320–500 mm. The proximal 22 mm of nails of diameters 10 and 11 mm are strengthened to diameter 11.5 mm. Diameters 12 mm to 13 mm. Diameters 13 mm to 14 mm. Diameters 15 mm the widening is not done.
The Femoral nails are designed in a way that the strengthening of proximal part does not have to be considered and the pre-drilling is according to the nail diameter (we generally pre-drill the diameter 1 mm larger than the chosen diameter of the nail). The nail is slightly bent so that it makes a circle stroke with a semi-diameter of 2.3 m. The part that is not strengthened has 3 mouldings. The nail is universal, can be used in the left and also in the right femur.

2) MEDIN Femoral intramedullary retrograde nails are hollow, have a circular section.
Distal part of 35 mm of nail of 10.5 mm is strengthened to diameter 11.5 mm. Nail of size 12.5 mm has the same diameter along the whole length. The nail is universal, the same for left and the right limb. The nail is hollow and that enables its implantation on a guide rod.
The nails sizes are: 200; 250; 300 and 350 mm in diameters of 10.5 mm and 12.5 mm.

3) MEDIN Tibial intramedullary nails are hollow without groove on the surface, have a circular section. MEDIN has at its disposal tibia nails of diameter 8 mm, sizes 240–360 mm, after 15 mm. They can be used as non-predrilled. Nails of diameters 9–12 mm, sizes 270–390 mm are also available. Proximal 75 mm are curved at the angle 13 ° in all sizes and lengths. Part of proximal 22 mm of nails of diameters 8; 9 and 10 mm is strengthened to diameter 11.5 mm. Tibial nails of 12 mm diameters are strengthened to diameter 13 mm in proximal part of 22 mm. The pre-drilling is according to the nail diameter, the strengthen proximal part does not have to be considered (is constructionally solved). E.g. we pre-drill about 1 mm more for a nail of diameter 11 mm, so the diameter of pre-drilled hole is 12 mm. The nail is universal, can be used in the left and also in the right tibia. A new type of an adjusting screw of 3.5 mm diameter is prepared for a tibia nail of 8 mm diameter. Because the opening of 5 mm diameter would make this nail too weak. A set of sizes relevant to chosen adjusting screws has to be available to customers.

Standard adjusting screws of 5 mm diameter are prepared for all types of MEDIN's intramedullary nails. They are supplied in sizes from 25 to 90 mm after 5 mm in three variations. First variation has the length of screw thread 20 mm on all sizes (the screw type is used in this instruction). The second variation has the screw thread up to the screwhead. A drilling sleeve nr. 4 is used for its implantation and a drill with diameter 3.5 mm. The third variation has a lowered thread profile, the thread is up to the screwhead. A drilling sleeve nr. 5 and a drill of diameter 4.4 mm is used. For only tibial intramedullary nails of 8 mm diameter, it is essential to use a screw of 3.5 mm diameter, drilling sleeve nr. 6, a drill with diameter 2.9 mm and sleeve nr. 7 for screw introduction – see MEDIN tibial intramedullary nails.

Colours of drilling sleeves and relevant drills
Sleeve 3 (129 69 1200) and drill ø 5 (129 79 4980) is distinguished by **white colour**.
Sleeve 4 (129 69 1210) and drill ø 3.5 (129 79 4990) is distinguished by **yellow colour**.
Sleeve 5 (129 79 8460) and drill ø 4.4 (129 79 8430) is distinguished by **red colour**.
Sleeve 6 (129 79 8450) and drill ø 2.9 (129 79 8440) is distinguished by **green colour**.

Set up of implant system
MEDIN's implant system of intramedullary nailing consists of one size of a nail, locking screws in proximal or distal part. It is recommended to close the nail with a stopper. It is possible to use a compression screw on tibial or femoral nails instead of a stopper and to perform static compression of fracture directly on the operating table.

RECOMMENDED PROCEDURE OF IMPLANT USAGE:
MEDIN FEMORAL INTRAMEDULLARY NAILS
At the research facility of the designer they prefer the patient to be put on an extension table on his/her back, the healthy limb is abducted in neutral rotation and the leg support is lowered to its maximum (hip dorsiflexion) while the injured leg is slightly abducted in neutral rotation, following the body's axis. Of course, the patient can be positioned on the side of his/her healthy leg, too. In cases of early operations, skeletal traction is not used. The position must allow X-ray control in the whole extent of the operated femur, place of the fracture in both projections, the rest of the femur in AP projection at least and place of distal locking in side projection (better in both).

At least approximate reposition of the fracture is performed and the operation area and the patient are prepared. The incision shall be led so as to extend the axis of the thigh, proximally from the major trochanter. The medullary cavity of the femur is opened by the perforator, dorsally at the medial edge of the major trochanter. If the full length of the cutting edge of the perforator is applied, an opening of 12.5 mm in diameter is created, which is needed for making a hole necessary for nailing without pre-drilling. By means of an awl of 5 mm in diameter, a channel is made in the spongiose bone of the proximal part of the femur and then an introducing wire with olive, 1150 mm long, is introduced through the fracture site, with additional reposition of the fracture if it is necessary.
The length of the nail shall be read by means of a radiolucent ruler or we can use the inserted introducing wire for this (verify the position of its distal end, measure the part extruding from the bone on another introducing wire of the same length and, finally, measure the remaining length).

If we use a nail without pre-drilling, a chosen nail of 10 or 11 mm in diameter is introduced and then the introducing wire is removed and the below described method of interlocking is applied. The pre-drilling must be performed gradually by means of the whole set of reamers, up to the selected size. It shall always be initiated with the reamer of 8 mm in diameter which is the only one designed for frontal reaming. The reaming shall be finished by a reamer of a diameter 1 mm bigger than is a diameter of selected nail. After irrigation (the introducing wire is replaced by a tube enabling irrigation and then re-introduced), the selected nail is assembled with a stirrup 1 (the same aimr is applied for the tibial and femoral nail, a different one for the reconstruction nail). The hammering can be performed directly on the specially modified part of the aimr 1 – the impactor. The method of screwing the impactor into the aimr is the same as that is applied with the hammer. It is recommended, though, to use (screw in) a hammer. In order to make the manipulation easier during the introduction of the femoral nail, the positioning of the guide rod of the hammer on the aimr is eccentric to the nail axis. The nail is introduced by strokes of the weight. Undesirable rotation of the nail during the introduction needs to be prevented. The introducing wire is then removed.

The proximal locking brings no difficulties thanks to using the aimr 1. It is reasonable to remove the hammer before applying the interlocking. You can introduce 1 or possibly 2 locking screws. Nail openings enable lateromedial introduction in the frontal plane only. The upper opening for the screw is oval-shaped and introduction of the screw by aimr enables both dynamic and static interlocking. Compression of the fracture site can be achieved right on the operation table by using a compression screw introduced to the upper end, axially into the nail cavity, with the last mentioned method of interlocking. The maximum extent of compression is 10 mm. (Distal locking must be performed and only one screw is introduced proximally into the oval-shaped opening of the nail. The compression screw shall lean on the proximal locking screw in the nail cavity and by further insertion it shall push the proximal fragment of the bone distally against the nail.)
In cases of more significant distraction of the fracture site, which might be caused by the introduction of the nail, it is better to bring fragments closer by an alternative method, keeping the aimr. The handle and weight of the hammer can be temporarily removed. Then we perform distal locking and reassemble the hammer. The approximation of the fracture site shall be achieved by reverse strokes of the hammer (like during extraction).

Via chosen opening of the stirrup 1, we introduce drilling sleeve assembled with a trocar for easier introduction. In the place of contact of the trocar tip we perform an incision of the skin of appropriate length, approximately 15 mm. We push the sleeve to the bone by hand first and then by gentle strikes of the hammer. By stroking the trocar 1 we create a depression in the bone so that the drill can be easily inserted. Trocar 1 is removed and the sleeve pushed in a way so that the toothed inner edge is fixed well on the bone. Drill sleeves 3 and 4 shall be inserted in sleeve. The sleeve 3 is for a drill of 5 mm diameter, the sleeve 4 is for a drill of 3.5 mm diameter. Now, we drill an opening through both the cortices with the 3.5 mm drill. We remove the inner sleeve 4 and drill the closer cortex by using the 5 mm drill through inner sleeve 3. We remove the inner sleeve and by using the depth gauge, we measure the needed length of the locking screw to be inserted. We can avoid any possible doubts concerning the position of the scale mark, against which we read the value, by setting the depth gauge to its zero position (completely closed), then the position of the mark is absolutely clear. The necessary length of the locking screw is read directly as the correction of the distance from the sleeve 2 to the outer cortex is done on the depth gauge!

Locking screws are self-tapping and can be introduced in this phase. Besides, tap 1 held in hand chuck can be applied as well. The same procedure can be used for the introduction of the second locking screw. The screw is introduced by a screwdriver, through sleeve 2.

Locking screws are not threaded all along and therefore, their extraction can be performed merely by using a screw extractor assembled with a screwdriver. If a screwdriver is equipped with a screw extractor, it enables screwing out at the same time.
Besides, drill sleeve 1 is supplied, which can be inserted into sleeve 2 and enables introduction of a K-wire of 2.5 mm in diameter in order to verify the position of the locking screw prior to commencement of the drilling process. Using the same

method, drill sleeve 1 can be assembled with the aiming device for distal locking (to be described further).
Prior to distal locking we remove aimr by loosening the screw of insertion handle with the socket wrench contained in the set of instruments and by subsequent extension of the body of aimr. If we don't intend to perform compression of the fracture right on the operation table with a compression screw, the nail cavity can be closed with a stopper.

The distal locking is more complicated. On the basis of preceding experience, locking "out of hand", with the help of an aiming device, has been selected. At first, the C arm of the mobile X-ray machine is set right into the axis.

The screen must show openings for locking screws exactly circular. Upon the incision of the skin, an aiming device with a trocar is introduced to the bone in the same way as in the proximal locking. The trocar having been removed, the aiming device is tapped to the bone. The position is checked by X-ray monitoring.

Sleeve 1 can be inserted into the aiming device now and the position is verified by introduction of 2.5 mm K-wire. Then sleeves 3 and 4 (eventually 5) are inserted and both the cortices are drilled with 3.5 mm drill. Sleeve 4 is removed and the closer cortex drilled with 5 mm drill. Sleeve 3 is removed. The length of the locking screw is measured with a depth gauge now. Tap 1 may be applied, if this is considered appropriate. The screw, having been selected, is introduced by means of a screwdriver now, still being led through the aiming device.

We shall introduce 1–2 distal locking screws, nail openings enable the lateromedial introduction in the frontal plane only.

A double aiming device is supplied for the individual choice of an surgeon. The double aiming device is assembled from the handle with the aiming device, sleeve and sleeves 1; 2; 3; 4; 5 and 6. The handle is loosely mounted on the aiming device, sleeve 4 is inserted into the protection sleeve of the aiming device and through it, a 2.5 mm K-wire, 160 mm long, is introduced through the distal locking opening of nail. Sleeve, sleeve 2 and sleeve 1 is then inserted into the second opening of the handle and, when they have been positioned correctly, the arresting screws of the handle are tightened. A K-wire can be inserted through sleeve 4 in order to verify the position. By means of sleeves 3–5, using the above-described method, an opening is prepared and the locking screw is introduced into the proximal one of the two distal nail openings. Upon extraction of the K-wire from the distal nail opening, the distal locking screw is introduced in the same way.

After irrigation, the incisions are sutured. At the site where the nail was introduced into the femur, a sucking drain is introduced on a regular basis; it shouldn't suck right from the medullary cavity. The incision is covered with a soft dressing and X-ray documentation is made. If skeletal traction has been used, it is removed.

RECOMMENDED TECHNIQUE FOR FEMORAL NAIL EXTRACTION
First, it is necessary to remove locking screws (if it has not been done during the healing process). The extraction is possible only by means of a screwdriver equipped with a screw extractor, which enables unscrewing and traction at the same time. The incision having been made, we put the screw extractor with a trocar to the screw to be extracted. The tip of the trocar should aim at the head of the screw. We remove the trocar and replace it by the screwdriver, which is to be assembled with the screw extractor by using pressure till it "snaps in".

The stopper or compression screw must be removed from the proximal opening. The extraction rod is screwed in to the inner thread of the proximal part of the nail. The nail is then driven out by striking the weight against the handle.

MEDIN FEMORAL INTRAMEDULLARY RETROGRADE NAILS
The placing of the patient is the same as placing during synthesis of common femoral nail.
The position on extension table or on a common operating table must enable to see not only the point of the fracture and the place of distal locking but also view in two levels. The knee flexion is 40–50 °.

- At a fracture of type C, first we make a reposition and from single incisions, trascutaneously, we stabilise the diacondylal fracture of femur by using of cancellous or cannulated screws.
- We make an incision in length of 5 cm across the ligamentum patellae. Trepanation of medullary cavity is made in the ceiling of intercondylal foss, in front of attachment of cross ligaments. We introduce flexible introducing wire and according to principles of intramedullary osteosynthesis we pre-drill a medullary cavity. The procedure is similar to pre-drilling of other types of nails. The X-ray check-up is necessary in two projections.
- Chosen nail is installed on the aiming device and it is introduced in a common way to the bone cavity. It is possible to check the position of indwelling by X-ray by means of grooves in the hub of the aiming device.
- Distal locking is a very easy matter. Sleeve 2 is put in appropriate openings of the aiming device and then sleeves are put too, according to a type of used locking screw:
 - Locking screw (partial thread) sleeve 4 is put in sleeve 2 and an opening of a drill of diameter 3.5 mm is pre-drilled through both of corticalis. Sleeve 4 is replaced with sleeve 3 and closer corticalis is drilled off by a drill of a diameter 5 mm.
 - Locking screw (thread up to the screw head) sleeve 4 is put in sleeve 2 and an opening of diameter 3.2 mm is pre-drilled through both of corticalis.
 - Locking screw 5 – Sleeve 3 is put in sleeve 2 and an opening of diameter 5 mm is pre-drilled through both of corticalis. The length of locking screw is measured by a depth gauge that is introduced through the sleeve 2. The locking screw is introduced through this sleeve with using of an appropriate screwdriver.
- Proximal locking is made by using the same above mentioned locking screws, by hand. It is made lateromedially when using shorter nails and ventrodorsally when using long nails.
- The nail is closed with a stopper. It is possible to fix the stopper in a way that the screw of the stirrup is removed from the aiming device, the aiming device is placed on the nail and the stopper is inserted through the hub of the aiming device. After irrigation, the incisions are sutured. A sucking drain is introduced regularly to the site where the nail was introduced into the femur; it shouldn't suck right from the medullary cavity. The incision is covered with a soft dressing and X-ray documentation is made.

Notes:
– Prior to nail insertion, when the nail and the aiming device is set up, it is necessary to check whether all openings to which the screws will be introduced match with the openings in the aiming device.
– Use a hammer during the nail insertion only adequately. The nail should normally be inserted only by hand as far as possible.
– Some implants are made of steel or titanium alloy. They can never be combined at one patient.

RECOMMENDED TECHNIQUE FOR RETROGRADE NAIL EXTRACTION
1. The stopper is removed. An extraction rod is led into the inner cavity of distal nail end. Weight and a handle are fixed on the rod. After that, it is essential to hold these parts because of possible overburdening of healed bones with moment of flexion. Caution! THE HEALED BONE COULD GET DAMAGED AGAIN!
2. The nail must be locked during placing the extraction rod due to prevention of its insertion into the bone cavity.
3. Remove last locking screws.
4. Take out the nail by strikes of the weight against the handle.

MEDIN TIBIAL INTRAMEDULLARY NAILS
The patient is placed on the extension table; traction is carried out more frequently skeletally, by drilling the calcaneum. Soft traction through a figure-of-eight shaped bandage on the anklebone is less common. At the research facility of the designer they prefer the patient to be placed on his/her back, with the operated leg following the body's axis, in neutral rotation, with the knee flexed so that the shinbone is positioned at a slight angle, distally.
The position must allow monitoring by the fluoroscopy intensifier, at least as follows: at the fracture site in both the projections, at the site of the proximal locking the AP projection is sufficient, the distal view depends on the type of the interlocking having been chosen but usually, the side projection is sufficient; the possibility of double-plane imaging being advantageous even in cases where locking screws are inserted only in one plane.

Under double-plane control by fluoroscope intensifier, at least an approximate repositioned of the fracture is performed. The operation area and the patient having been prepared, a longitudinal incision is led above the centre of the ligamentum proprium patellae in its full length. The ligament is halved, both the halves are drawn back and held in place by raspatories being driven into the bone. By means of the perforator contained in the set of instruments, the medullary cavity is then opened above the upper edge of the tuberositas tibiae. (When buying the set of instruments, you can choose one of two available types of perforators; the case contains a place for one of them.) If the perforator is driven at the full length of its cutting edge, the opening (in both the types used) has a diameter of 12.5 mm. Then the introducing wire is inserted (in the case of a tibia 950 mm long, with an olive when pre-drilling is applied or, in nailing without pre-drilling an introducing wire is driven in, possibly after the 5 mm awl has been inserted) into the distal fragment with additional reduction of splinters being necessary.

Sometimes it is more reasonable to localize the opening for trepanation slightly medially from the central line. Then the skin incision can be done at the medial edge of the ligamentum proprium patellae and the ligament drawn back laterally.
If a 12.5 mm opening has been made by means of the perforator contained in the set of instruments, it is not necessary to respect the strengthened proximal end of the nail and the pre-drilling is done according to the diameter of the chosen nail. In a case of nailing without pre-drilling, the chosen intramedullary nail is driven directly. We can use diameters 8; 9 or 10 mm.
During pre-drilling, we use the whole set of cutters up to a diameter 1 mm bigger than the diameter of the chosen nail, in 0.5 mm steps. The cutting should always start with the 8 mm cutter as this cutter is the only one designed for frontal cutting.

The nail's length is read by means of a radiolucent ruler or we can use the introducing wire inserted for this (we verify the position of its distal end and measure the remaining part that is protruding from the bone at another introducing wire of the same length and read the remaining length).

After irrigation (the introducing wire is replaced by a tube enabling irrigation and then re-introduced), the selected nail is assembled with aimr with the help of a screw (the same aimr is applied both for the tibial and femoral nail, a different one for the reconstruction nail). The hammering can be performed directly on a specially modified part of the aimr – the "impactor". The method of screwing the impactor into the aimr is the same as the one applied with the hammer. It is recommended, though, to use (screw in) a hammer. The guide rods of both the hammer and impactor should be tightened well when used so that they do not get flexurally stressed. In order to facilitate manipulation during the introduction of the tibial nail, the positioning of the guide rod of the hammer on the aimr is eccentric to the nail axis. The nail is then driven in with strikes of the weight. Undesirable rotation of the nail during the introduction needs to be prevented. The introducing wire is then removed.

The proximal locking causes no difficulties when using the aimr. It is reasonable to remove the hammer before applying the interlocking. You can introduce one or possibly two locking screws. Nail openings allow mediolateral or better, latero-

medial introduction only. The upper opening for the screw is oval-shaped and introduction of the screw through the aimr enables both dynamic and static interlocking. Compression of the fracture site can be achieved right on the operation table by using a compression screw 1 introduced at the upper end, axially into the nail cavity, with the last-mentioned method of interlocking. The maximum extent of compression is 10 mm. (Distal locking must be performed and only one screw introduced proximally, into the oval-shaped opening of the nail. The compression screw shall rest on the proximal locking screw in the nail cavity and, by further insertion, it shall push the distal fragment of the bone to the proximal one.)

In the process of inserting the nail, the fracture site may be distracted and the fragments may need to be brought closer by more than 10 mm. In such cases it is reasonable to leave aimr at its place and, possibly to temporarily remove the handle and weight of the hammer. Then the distal locking can be performed and the hammer reassembled. The approximation and certain compression of the fracture site shall be then achieved by reverse strokes of the hammer (like in extraction). This can be combined with the technique described in the preceding paragraph.

Introduce the drilling sleeve 2 assembled with a trocar via a chosen opening of the aimr for easier introduction. At the contact place of the trocar tip we perform an incision of the skin of appropriate length, approximately 15 mm. We push the sleeve to the bone by hand first and then by gentle tapping the hammer. By stroking the trocar we create a dent in the bone so that the drill can be easily inserted. Trocar shall be removed and the sleeve pushed in a way that the toothed inner edge is mounted well on the bone. Drilling sleeves 3.4 and eventually 5 shall be inserted in sleeve 2. Sleeve 3 is for a drill of 5 mm in diameter. Sleeve 4 is for a drill of 3.5 mm in diameter. Now, we drill an opening through both the cortices with the 3.5 mm drill. We remove the inner sleeve 4 and, using the 5 mm drill, drill the closer cortex. We remove sleeve 3 and by using the depth gauge, we measure the needed length of the locking screw to be inserted. We can avoid any possible doubts concerning the position of the scale mark, against which we read the value, by setting the meter to its zero position (completely closed), then the position of the mark is absolutely clear. The necessary length of the locking screw is read directly as the correction of the distance from the sleeve 2 to the outer cortex is done on the meter! Locking screws are self-tapping and can be introduced in this phase. Besides, tap 1 held in hand chuck can be applied as well. The same procedure can be used for the introduction of the second locking screw. The screw is introduced by a screwdriver, through sleeve 2. Locking screws that are not threaded all along, therefore their extraction can be performed simply by using a screw extractor assembled with a screwdriver. The set containing a screwdriver and a screw extractor enables unscrewing a screw and use traction at the same time.

Besides, drill sleeve 1 for nail interlocking is supplied in the set of instruments. It can be inserted into sleeve 2 and used for introduction of a K-wire of 2.5 mm in diameter (160 or 300 mm long) in order to verify the position of the locking screw prior to commencement of the drilling process. Using the same method, drill sleeve 1 can be assembled with the aiming device for distal locking (to be described further).

Prior to distal locking we remove aimr (in case we do not intend to bring closer the fragments by means of the technique described above), by loosening the screw of the aimr with the socket wrench contained in the set of instruments and by subsequent extension of the body of aimr. We can perform compression of the fracture right on the operation table with a compression screw or close the nail cavity with a stopper.
The distal locking is more complicated. On the basis of preceding experience, interlocking "out of hand", with the help of an aiming device, has been selected. At first, the C arm of the mobile X-ray machine is set right in the axis. The screen must show openings for locking screws exactly circular. Upon the skin incision, an aiming device with a trocar is introduced to the bone in the same way as in the proximal locking. The trocar having been removed, the aiming device is tapped to the bone. The position is checked by X-ray monitoring.

Sleeve 4 can be inserted into the aiming device now and the position verified by introduction of a 2.5 mm K-wire. Then sleeves 3; 4 and possibly 5 are inserted and both the cortices are drilled with a 3.5 mm drill. Sleeve 4 is removed and the closer cortex is drilled with a 5 mm drill. Sleeve 3 is removed. The length of the locking screw is measured by a depth gauge now. Tap 1 may be applied, if this considered appropriate. The screw, having been selected, is now driven in by means of a screwdriver, still being led through the aiming device.

We shall introduce 1–3, eventually 4 distal locking screws, with nail openings allowing for double-plane introduction.

A double aiming device is supplied, to be used at the discretion of the surgeon. The double aiming device is assembled from the handle and sleeves 1; 2; 3; 5 and 6. The handle is loosely mounted on the aiming device, sleeve 1 is inserted into the protection sleeve of the aiming device and through it, a 2.5 mm K-wire, 160 mm long, is introduced through the distal locking opening of a nail. Sleeves 1 and 2 are then inserted into the second opening of the handle and, when they have been positioned correctly, the arresting screws of the handle are tightened. A K-wire can be inserted through sleeve 1 in order to verify the position. By replacing sleeve 1 with sleeves 3 and 4, eventually 5 and 6 an opening is prepared and the locking screw is introduced into the proximal one of the two distal nail openings. Upon extraction of the K-wire from the distal nail opening, the distal locking screw is introduced in the same way.

The incision above the ligamentum proprium patellae is irrigated and all the incisions sutured.

X-ray documentation usually contains images of selected views from the intensifier screen, however, with the necessity for a classical X-ray double-plane picture made on a large-scale format for precise determination of the axis of the synthesized bone. If skeletal traction has been used, it is removed now. The operated leg shall not be fixed, upon placement of a soft dressing it is reasonable to put an appropriate elastic bandage above the knee.

RECOMMENDED TECHNIQUE FOR TIBIAL NAIL EXTRACTION
First, it is necessary to remove locking screws (if it has not been done during the healing process). The extraction is possible only by means of a screwdriver equipped with a screw extractor, which enables simultaneous unscrewing and traction. Having done the skin incision, we introduce the screw extractor with a trocar to the screw to be extracted. The tip of the trocar should aim at the head of the screw. We remove the trocar and replace it with the screwdriver, which is to be assembled with the screw extractor by using pressure until it "snaps in".

The stopper or compression screw must be removed from the proximal opening. The rod hammer is screwed in to the inner thread of the proximal part of the nail. The nail is then driven out by striking the weight against the handle.

RECOMMENDED TECHNIQUE OF OSTEOSYNTHESIS WITH INTRAMEDULLARY NON-PREDRILLED NAILS.
Certain sizes of the line of the MEDIN intramedullary nails are designed for insertion without pre-drilling in tibial nails - those diameters are 8; 9 and 10 mm and for femoral nails diameters of 10 and 11 mm. Reaming of the intramedullary cavity does not take place here, of course, but otherwise the technique is in no way different from the above-mentioned methods of osteosynthesis.

MEDIN REAMERS OF INTRAMEDULLARY CAVITY
The line of MEDIN reamers serves for pre-drilling of the tibial medullary cavity. The complete line starts with a reamer of 8 mm in diameter and continues up to 16.5 mm in 0.5 mm steps. Up to 16 mm they belong in the standard set, the 16.5 mm size or, alternatively, a rigid 16.5 mm drill is to be purchased individually. A quick coupler enables connecting of both the reamers and drill to commonly used driving units.

It is not recommended to use the drill in reverse! "Unreeling" of the flexible spindle of the reamers might happen even during slightest weighting.

We request using the reamers only with original conductors MEDIN, diameter 3 mm and length 950 and 1150 mm, ended with an olive. It is prohibited to use the reamers without these conductors. Reaming operation has to start with the cutter of 8 mm, this is the only one specially adjusted for frontal cutting. It is not recommended to omit some sizes of cutters; consequently it is necessary to do the reaming after 0.5 mm.

To introduce reconstruction nail it is essential to make wider proximal enter opening and a canal in a medullary cavity to the diameter at least 15.5 mm in the distance of 80 mm (when the cutting of reamer with the whole set is finished). The last size of reamers is used for it. The manufacturer is able to supply the reamers with diameters: 15.5 mm; 16 mm and 16.5 mm.

RECOMMENDED METHOD OF CLEANING
1) Mechanical cleaning by water and a brush.
2) Rinsing instruments by compressed water, the best is to install a nozzle on a water valve and by means of this pressure to wash the instruments.
3) To immerse the instruments into disinfection solution at least for 10 minutes. We recommend 20–30 minutes. A recommended disinfectant is KORSOLEX plus.
4) Repeated rinsing by the stream of compressed water until the water is clean. This is possible in conventionally supplied, compressive washing machines when the instruments are connected to the nozzles.
5) Cleaning the instruments by pressed air.

NOTE:
While any using and even during cleaning the flexible spindle of the reamer must not be bent too much in particular bending in small diameter. Allowed bending radius is given by bending of reamer which is put on original conductor specified by the producer for the pre-drilling during which excessive weighting is not allowed. When the work is finished all the dirt and soiling should not get dry on the reamers. The cleaning would then need more attention and there would be no guarantee that the instruments would be cleaned to perfection.

INSTRUMENTATION
The instruments are made of austenitic or martensite stainless steel that is used for these purposes for long time. They are supplied non-sterile but clean. The implants from the instrumentation must be sterilised before use. The instrumentation is also supplied non-sterile but clean. Recommended sterilisation is in steam steriliser under 121 °C, overpressure 105 kPa for at least 20 minutes. The implants and also the instruments are packed separately in PE bags that have self-sticking labels on them.
The instruments of the instrumentation can be cleaned repeatedly, disinfected and also sterilised, cared according to PL0138 – instruction of use. The implants must be cared according to general enclosed leaflet for implants PLO088. The implants are disposable, for one patient only and for one stabilisation of damaged bone, repeated use is prohibited. This fact is mentioned in enclosed leaflet concerning all implants. Single instruments can be used repeatedly for the time set by a doctor (if not mentioned differently).
The instruments and implants must be liquitated in hospitals in the same way as other infected material, see PLO088.