**CZ**

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### VAROVÁNÍ

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**EN**

Information on use PL0148  
**READ CAREFULLY!**

## MESH STRENGTHENER

### WARNING

It is possible to use gold-plated mesh for the patients who are allergic to nickel.

The dentures that are not fixed but can be taken out of the mouth must not be mechanically damaged, otherwise that could cause swelling.

### Restriction of repeated usage:

The products are disposable. The products are determined for one introduction into the prosthesis only and the manufacturer does not recommend repeated using. Other using can cause functionality loss and it can also cause contamination or illness.

### INSTRUCTIONS

#### Where to use it:

The meshes are anatomically shaped plates, made of steel and are intended for stabilisation of upper complete or partial replacements. The holes in the mesh enable a fix connection of the resin layer over and under the mesh. The mesh increases the strength and lifetime of the prosthesis.

Cut the mesh to the required size (attention to the changes during consequential shaping – shortening) and pre-shape on the readymade duplicate of the model. The mesh is perfectly shaped according to the model in the cell press. Then cut off the excess parts. Further proceed in common way. During prosthesis making, pay attention to the mesh, so that it would be surrounded

by plastic from all the sides. The mesh border of the complete prosthesis should reach approximately to the half of the height of the alveolar ridge. At the partial replacements, it should reach further than to the half of the ridge, but it is necessary to keep sufficient distance from the interdental papillas of remaining teeth. Dorsally, the mesh border should end ca 2 mm before the back edge of the replacement.

### Keeping and transport:

Keep the meshes in solid packing to prevent from damaging.

### Storing:

Store the packed products in a dry and clean environment without extreme temperatures and humidity, out of the reach of direct sunlight.

### Disposal of waste:

The patient disposes the dentures on his/her own, or the dentist (it ends up in mixed waste).

### Other information:

Company MEDIN, a.s. validated above mentioned procedures. The procedures are qualified for preparation of medical devices of repeated use.

The user is responsible for ensuring of necessary materials, equipment and employee training in compliance with these set procedures so that required result could be reached.

The user can also apply other procedures but the effectiveness of these alternative procedures must be verified by a qualified laboratory.

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