MOMENTOVÉ INSTRUMENTY

ÚČEL POUŽITÍ

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TORQUE TOOLS

Intended use

The torque tools are used during the implementation of angle stable bone screws (hereinafter referred to as the screws) which must be tightened by means of defined torque. The use of the torque tools prevents the tightening of screws into the splint by means of torque which could cause cold weld between the thread surfaces contact areas, and significantly complicate the implant extraction. The tools are preset by the manufacturer to a specific torque value. The torque tool indicates a correct torque by means clicking during the screw tightening. The tool can indicate the correct tightening sooner, such as during the screw implementation through the splint into high quality bone material which requires higher torque for the screw tightening than defined by the instrument. In such case use a standard screwdriver, and use the torque tool for tightening the screw into the splint. The torque tools can be combined with replaceable ends and a holder.

Review of tools:

| Tool | Order number | Application | Torque preset |
|--------------------|--------------|--|---------------|
| Torque screwdriver | 129 69 2810 | Angle stable screws 3.5 and 3.5/2.7 mm | 1.5 Nm |
| Torque limiter | 129 69 5126 | Angle stable screws 3.5 and 3.5/2.7 mm | 1.5 Nm |
| Torque limiter | 129 69 5121 | Angle stable screws 5 and 5/3.5 mm | 4 Nm |
| Torque limiter | 129 99 1175 | Knee replacement | 2,75 Nm |

Note: The principal, procedure, and potential problems of angle stable splints are specified in the attached leaflet PL0180.

WARNINGS

Do not use the tools together with any electrical and pneumatic drills, although the construction allows it. Avoid agents containing chloride ions during cleaning and disinfection; they can cause corrosion.

Restriction of reoccurring processing

Reoccurring processing has a minimum impact on the tools. The end of service life is commonly specified by wear and tear and damage during use.

INSTRUCTIONS

Application sites:

Specialist workplaces of healthcare providers, particularly traumatology, surgery, and orthopaedics.

Storage and transport:

Protect the tool from deformation and scratching; avoid contact with chemical agents and other metal materials during handling, except for those specified in the intended use.

Preparation for cleaning:

No special requirements.

Cleaning – automatic:

The tool must be inserted into the device for washing and disinfection. The washing cycle and rinsing

must be completed. During removal, check if visible contamination was removed. If required, repeat the whole cycle or clean the tool manually.

Cleaning – manual:

Excessive contamination of the tool requires rinsing under warm water. Use the brush with detergent solution to clean the whole surface. Repeat thorough rinsing under running water and leave to dry.

Disinfection: We recommend using device for washing and thermal or thermal-chemical disinfection, or manual washing of the tools after chemical disinfection by means of agents with virucidal effects; both stages can be combined by means of disinfection agents with combined washing effects. The disinfection agent can be used in compliance with the instructions specified on the label. Adhere to the procedure specified by the manufacturer for the preparation of disinfection and washing solutions.

Recommended agents: Sekusept pulver, Neodisher SeptoMED

Drying:

If drying is a part of the tool cycle for washing and disinfection, temperature should not exceed 120 °C.

Maintenance:

Damaged tool must be rejected.

Inspection and function testing:

Torque inspection is performed once in twelve months. Unchecked tool must not be used. Check the tool visually for damage and wear and tear. Damaged and worn tools must be rejected.

Packaging:

It is possible to use a standard packaging material. Ensure the package is sufficiently large in order to place the tool inside avoiding problems with the packaging lids. The tool must be placed in the package or sterile containers which are then packed.

Sterilization:

We recommend to sterilize the tools separately from tools made from other materials (titan, plastic materials, ceramics, etc.). We recommend sterilization by means of humid heat in steam sterilizer equipped with antibacterial filter at temperatures 121 °C, overpressure 205 kPa for 20 minutes, or at temperature 134 °C, overpressure 304 kPa for 10 minutes.

Hot air sterilization is completed in devices with enforced air circulation at temperature 160 °C for 60, or at temperature 170 °C for 30 minutes, or at temperature 180 °C for 20 minutes, unless specified otherwise. Hot air sterilizer opens upon the cycle completion after cooling to at least 80 °C.

Make sure the maximum sterilizer load capacity is not exceeded during the tool sterilization in one autoclave cycle.

Storage:

Use a sieve for storing, transporting, and sterilizing the tool.

Other information:

The joint-stock company MEDIN validated the above mentioned methods of cleaning, disinfection, drying and sterilization. Through the validation were all procedures recognized as qualified for the preparation of medical instruments for repeated using. The user is responsible for ensuring of necessary materials, equipment and staff training in accordance with these specified procedures to reach the required effect. The user can use some other methods, but the efficiency of these alternative methods must be verified by specialized laboratory.

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