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**ANGULARLY STABLE PLATES WITH SCREWS 5 mm**

Caution: Only a specialist can work with the implants.

**IMPLANTS FOR FRACTURES WITH SMALL FRAGMENTS**
**LIST OF PRODUCTS**
**Angularly stable plates for screws with a shank of a diameter 3,5 mm and 5 mm (head Ø 7 mm)**
**Straight plate L × X holes**

Implant steel: 397 129 77 9330-9410  
Titanium: 397 129 77 9333-9413

**Locking screws**
**Locking screw Ø 5/3,5 mm; self-tapping**

Implant steel: 397 129 77 8231-8451  
Pre-drilling with a drill Ø 2,9 mm

**Locking screw Ø 5/3,5 mm; self-tapping**

Titanium: 397 129 77 8234-8454  
Pre-drilling with a drill Ø 2,9 mm

**Locking screw Ø 5 mm; self-tapping**

Implant steel: 397 129 77 7981-8211  
Pre-drilling with a drill Ø 3,9 mm

**Locking screw Ø 5 mm; self-tapping**

Titanium: 397 129 77 7984-8214  
Pre-drilling with a drill Ø 3,9 mm

**Locking screw Ø 5 mm cannulated Ø 1,6; self-tapping**

Implant steel: 397 129 77 8681-8901  
Pre-drilling with a drill Ø 3,9 mm

**Locking screw Ø 5 mm cannulated Ø 1,6; self-tapping**

Titanium: 397 129 77 8684-8904  
Pre-drilling with a drill Ø 3,9 mm

**PRODUCT DESCRIPTION**


Angularly stable plate may have one of the three basic types of holes for plate fixing to the bone with the help of screws, namely in their optional combination. However, it always contains at least one locking hole.

**Locking hole**

Self-tapping or self-drilling screws can be introduced into the bone and "locked" by fastening of the conical thread on the screw head in the same screw in the plate hole. Such connection is towards the plate angularly stable, the whole system works on the principle of inner fixator. This solution helps to prevent from following problems:  


- primary loss of reposition
- secondary loss of reposition, in the first place in cases of fragmentation fractures without sufficient bone support, low quality or osteoporotic bone
- periosteum compression and following worsening of blood supply of cortical bone

Use locking screws in this hole. It is also possible to use standard cortical screw with ball head.

**Standard oval hole**

Standard oval hole enables to the surgeon to use compression screws for fastening of broken fragments and fracture fixation by axial compression. The screws can be introduced under different angles, lengthwise or transversally.

Do not insert locking screws in this hole!


**Autocompression locking hole**

Autocompression hole enables dynamic compression in lengthwise axis of the plate at eccentric screw placing. Holes are not regular, compression is always effective in one direction only, from that reason the holes are placed in a way so that the compression would work towards the plate centre. When using locking screws the hole behaves in the same way as locking hole.

**System of angularly stable plates**

The system is based on the principle of the inner fixator. By means of the tight connection "screw – plate", the angularly stable connection is achieved between the plate system and the bone. This transfers part of the weight from the screw-bone connection to the screw-plate connection. Thanks to that it is possible to use the system successfully even in a less quality and an osteoporotic bone.

If the angularly stable plate is used with locking screws, the plate and the screws make tight system together. The plate must not be pushed against the bone so that the bone would still get sufficient supply of blood. Locking screws can be introduced monocularly without enfeebled the strength of the configuration or the fracture fixation, because the screws are firmly anchored in the plate.

**PRINCIPLES OF FIXATION**
**1. Only compression screws**

It is about a commonly used method of fracture fixation with the help of plate and compression screws. Fragments of a bone are fastened to the plate, the plate is closely pushed against the bone and stability of the whole system is reached by bone pressure on the bottom part of the plate with the help of compression screws.

Shape of the plate is crucial for the correct function of osteosynthesis.

**2. Only locking screws**

If you use only locking screws then you make a system that works on a very similar principle as outer fixator.

When using only locking screws, the plate shape is not what matters but high-quality fracture reposition before tightening of the locking screws.

**3. Combination of locking and compression screws**

- Use the plate on the base of the common plating principle; first perform fracture reposition and fixation by compression screws. (You can use the spacer to keep the space between the bone and the plate). Then screw the locking screws to reach the angularly stable plate fixation towards the bone.
- After the fracture reposition and plate fixation by locking screws, it is possible to fasten the separated fragment (e.g. multi-fragmentation or splintery fractures) to the fixed part of the bone by compression screw.
- After fixation of one fracture fragment by locking screws, it is possible to push the second fragment to the first one with compression screws and eventually after that to fix even the second fragment with other locking screws.

**Note:** Do not use locking screws closely to both sides of the fracture because the strain, to which the plate is exposed, would be concentrated only to the space between these screws and the plate could get broken. It is generally recommended to leave from 1 to 3 holes empty above the fracture area.

**Note:** To use the locking plate system correctly, it is necessary to understand the principles of fixation and their biomechanical properties. Draw your attention to pre-operation planning.

**INDICATION**
**Straight plate**

Straight plate with holes is suitable for fracture fixation of straight bones, hollow bones and the collar bone. It is possible to shape the plate before using so that it would copy the anatomical shape of the bone better. The plate is strong enough so that it is also suitable for heavier patients and bone fractures where a higher load can be expected.

**CONTRAINDICATION**
**Straight plate**

See general information about implants using – Instructions nr. PL0088.

**Combination with other products**

It is possible to combine plates with screws and wires of MEDIN, a.s., providing the instructions and producer's recommendations are followed, particularly compatibility of the materials. (it is not allowed to combine products of different materials)

The products are not determined to be combined with products of other producers.

**PLATE INTRODUCTION AND LOCKING WITH STANDARD SCREWS**

In case of using only standard i.e. unlocking screws, the operating procedure is the same as during plating of standard plates for osteosynthesis.

**1. Reposition the fracture**

Reposition the fracture under X-ray control and fixate it with K-wires MEDIN or forceps (397 129 99 0040).

**2. Shape the plate**

Get to know the fracture place really well. For correct function, it is necessary to shape the plate anatomically correct in advance, mainly for internal-articular fractures.  
Use only allowed shaping instruments, those that are identical with instruments for standard plates.

**3. Insert guide sleeve locking (Alternative)**

Screw the guide sleeve locking carefully into the thread of chosen plate hole until it is tightly hold in the thread. Use the guide sleeves locking for manipulation and correct plate placing.

**4. Place the plate**

Place the plate on the bone and fix it temporarily with forceps (397 129 99 0040) or K-wires MEDIN.

**5. Pre-drill the hole for screw**

Pre-drill the hole for the screw with the help of standard instruments for plate osteosynthesis. Screw can be introduced under angle either lengthwise or transversely with the plate axis.

**6. Insert the screw**

Introduce chosen standard screw neutrally or in case of autocompression holes you can place the screw eccentrically and to reach dynamic compression this way.

**7. Make sure everything is all right**

Check that the screws are tight. Remove instruments necessary for temporary plate fixation.

**8. Plate removal**

The plate removal procedure is the same as procedure of standard plates. If the hexagon head of the screw gets stripped or during other problems with screws removing use the instrumentation for removal of MEDIN screws.

**PLATE INTRODUCTION AND LOCKING WITH LOCKING SCREWS**

Despite of the fact that locking screws are introduced bimorally in most cases, it is possible to introduce them even monocularly in quality and strong bones. Monocortical introduction of the screws is mainly used in case when the screw is not desired to go through the whole bone e.g. in the area of epiphysis.

Pay attention to locking screws, they do not behave as standard compression screws. It is necessary to realize that during fastening the screw and locking the screw's head in the plate hole, the screw is strongly fixed and the force compression of bone fragments does not happen. Therefore, careful anatomical bone reconstruction is necessary (mainly in intra-articular fractures) with the help of compression (standard) screws, K-wires or cannulated screws. If the intra-fragment compression is required, using of standard screws is also necessary.

**1. Reposition the fracture**

Reposition the fracture under X-ray control. If essential, fixate it with K-wires MEDIN or forceps (397 129 99 0040).

**2. Shape the plate**

Shape the plate before operation suitably. Use only allowed shaping instruments, those that are identical with instruments for standard plates.

It is not necessary to shape the plate exactly on the bone for correct and quality reposition. Quality and stability of the connection does not depend on the contact of the plate with bone surface. This is mainly suitable for mini-invasive operation methods, when it is not possible to shape the plate in advance.

When bending, prevent from deformation of the holes with thread, they could get damaged. Do not use locking screws or spacers to protect the thread during bending, they could get stuck in the hole, and thus their removal would be impossible or damage of the thread would happen. Do not use locking drilling sleeves for bending.

**3. Insert the spacer**


To minimize the contact of the plate with the bone, it is possible to place the spacer in the hole with thread before fixation of standard screws, it ensures the space between the plate and the bone 2 mm. It is possible to remove it again after plate securing with locking screws.

**4. Insert guide sleeve locking**

Screw the guide sleeve locking carefully into the thread of chosen plate hole until it is tightly hold. Use the guide sleeves locking for manipulation and correct plate placing.

**5. Place the plate and fix it temporarily**



Place the plate where it should be with the help of guide sleeves locking and fix it temporarily with the help of cortical screw in the oval hole, or possibly with K-wires MEDIN or forceps.

Correct position of the screws is necessary for right function, mainly screw alignment with the plate hole. To determine correct screw position you can use K-wires. Screw locking sleeve into required plate hole and insert reduction into it for wires and introduce the wire Ø 1,5 mm. Use RTG to picture the wires – they represent final position of locking screws.

**6. Pre-drill the hole**

Remove the reduction and leave the drilling sleeve with wires in the place. With the help of this sleeve, pre-drill the hole for screw of Ø 5 mm or cannulated screw of Ø 5 mm with the help of a cannulated drill of Ø 3,9 mm. If you decide to use screws of Ø 3,5 mm, put a reduction with inner Ø 3 mm into the drilling sleeve and pre-drill a hole for this screw with the help of a drill of Ø 2,9 mm. Drill only until first cortex then you can screw the screw easily thanks to the self-tapping thread on the tip. Carefully with a drill, pre-drill hole in the bone with the help of attached guide sleeve. Use the electric drill for mechanical pre-drilling. Use the hand-handle of the drill for more exact control above depth of drilling and to prevent the soft tissues from damage on the opposite side of the bone.

**7. Unscrew the guide sleeve**

Unscrew the guide sleeve carefully. In case that it is not possible to unscrew the sleeve only by hand from some reason, you can help yourself with guide sleeve for wires – insert it into the hole of the drilling sleeve handle and use it as an elevator. Use the wire hole for screw depth determination with gauge.

**8. Measure the hole depth and choose the correct screws**

Measure the depth of the pre-drilled hole with the help of the depth gauge (397 129 99 0040).

Perform it from the upper surface of the plate – the lengths of the screws are mentioned including the head. Choose the correct length of the screw carefully. If the bone is strong and of a good quality, it is possible to use the screws monocularly. For osteoporotic or less quality bone, we recommend using of bicortical screws. For proper system function it is the screw density important (proportion of the number of the screws to number of the holes in the plate). Ideal proportion is 0,4–0,5. In contrast with standard plates, using of screw in each plate hole is not recommended.

Nevertheless, plate fixation proximally or distally towards the fracture place is still very important.

At least two screws should be used monocularly for fixation of each main fragment, we recommend using 2 or 3 screws into main fragments, and in case fixation some of the screws were not optimal. To enhance the quality and stability of the system we recommend using at least one screw bicortically.

**9. Insert the screw**

Before introduction of the first locking screw, make sure that the plate is firmly fixed – when fastening the screw the rotation could happen – the plate round the screw axis. Make the screws tight with the help of hexagonal cannulated screwdriver.

Pre-drill the holes for locking screws with the help of the guide sleeve. For proper plate function, it is necessary to introduce the screws in the plate in a way that the lengthwise axis of the screw is identical to the axis of the corresponding hole in the plate. At a divergence smaller than 5° the connection of the screw's head and the plate is decreased to 70%.

**10. Make sure everything is all right**

Check that the locking screws are tight. Remove instruments necessary for temporary plate fixation.

**11. Plate removal**

When removing the plate, loose all the screws first, after that start to take them out one after the other. In this way, you prevent from possible rotation of the plate round the axis of the screw that was taken out as the last one. If the hexagon head of the screw gets stripped or during other problems with screws removing, use the instrumentation for removal of MEDIN screws.

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**CAUTION****PRECAUTION:**

1. The implant may be never repeatedly used, previous stress can cause defects that can reduce its lifetime.
2. The implant cannot be unlimitedly resistant to the level of functions which are in accordance with normal healthy bone. It can be used as long as the bone is not healed, the implant for fracture fixation at active younger patients is usually removed after this healing time. However at the older patient the surgeon can consider the implant leaving for fracture fixation after the healing than risk the additional surgical intervention of its evulsion even if the implant loses its function yet.
3. The surgeon must be acquainted with the operation process by using of specified instruments.
4. Use only the specified instruments for surgical procedure; it mustn't be injured to be used well.
5. Consider the use of implants for patients where you aren't sure that they will cooperate with doctor after operation (e.g. drunks, drug-addicts, mental sick).
6. Protect the implant from deformation and scratching, avoid the contact with other metal materials and chemicals by manipulation.
7. If the adjustment of implant shape is necessary, make it only with specified way.
8. Don't use various material-different implants for one patient, because the contact of two various metals could speed up the process of corrosion what could result in earlier drop of the implant effect and lifetime, further in increased volume of metal alloys that release in the patient's body and could cause a toxic reaction.
9. Decide the possible untimely implant removing if you surmise the bad bone healing or that the implant can lose, bend, break in the bone, because so it get to heavy patient's injury.
10. Dispose of all shapes and sizes of implants in the medical workplace which use can be supposed.
11. Take into consideration the temporary preventive means that reduce the weight loading for the purpose of minimization of early loading of the mean for fracture fixations as long as the evident bone junction isn't possible to see.
12. Decide the implant removing for internal fixation with taking account of factors that are risks for patients in consequence of additional surgical procedure and troubles during and after implant removing.
13. Patients which smoke or usually use product with nicotine have escalated number of uncured fractures. Patients like this must be warned of this fact and eventual consequences. Occurrence of implant failure is rather probable at heavy, physically active person.

**Information for patient:**

1. Warn the patient that he must exactly follow the instructions of nursing doctor.
2. Instruct the patient of operation process, possible risk and rehabilitation program.
3. Instruct the patient when he can fully load the bone with internal fixation in order to avoid the implant damage.
4. Warn the patient that the bone is predisposed to injury during treatment above all in points when the implant ends.
5. Announce to the patient when the implant is supposed to be removed.
6. Instruct the patient that he should abide adequate activity level after implant removing because it can get to repeated fracture before the bone is wholly healed and the bone cavity filled.

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**Possible side effects:**

1. Non-executed junction of fracture ends, osteoporosis, repressed revascularisation and bad ectostosis result in loosening, bending and breaking of implant or early loss of firm bone fixation.
2. Bad junction of fracture point following from wrong equalization.
3. Heightened reaction of fibrous tissue around the fracture point because of unstable comminuted fractures.
4. Early or late infection, both deep and surface.
5. Avascular necroseon.
6. Shortening of affected bone side (fracture)
7. Subclinical injury of nerves can appear as a result of surgical trauma.
8. Patient's metal sensitivity.
9. Risk of faulty diagnosis by dint of diagnostic devices CT or MRI where the occurrence of steel implant in the given locality causes wrong view.
10. The heavy complications can occur with every surgical operation for aftermath elimination of multiple injuries. These complications include except of that also urogenital, gastrointestinal and vascular injuries including thrombosis, bronchopulmonary injuries including embolism, heart failure and death.

**Limit of repeated treatment:** The implants are meant for disposable use.

Implants can be repeatedly disinfect and re-sterilized.

**INDICATION**

The implant is meant for holding of fractured bone surface in mutual contact for the purpose of fracture healing facilitation. The indications of concrete implants are mentioned on relevant separate package leaflet.

**CONTRAINDICATION**

1. Deficient quantity or quality of bone that can hinder the correct bone fixation.
2. Whichever fully developed or supposed latent infection.
3. Previous course of infections.
4. Whichever patient's mental disorder or neuromuscular disease that would lead to unacceptable risk in connection with fixation failure or with complications in the postoperative phase
5. Aggravated vascularity that would hinder the necessary blood distribution to the fracture or operating point.
6. Malignancy in the fracture point.
7. In all these cases when the instrument could come into conflict with anatomical or physiological structures.

**INSTRUCTIONS**

**Attention:** Before package opening it is necessary to check if it is not damaged. If it is damaged it is necessary to give implant back to manufacturer.

**Where to use it:** Professional workplaces of health care providers, mainly surgical and orthopaedic.

**Storage and transport:** Protect the implant from deformation and scratching; avoid the contact with other metal materials and chemicals during manipulation.

**Preparation for cleaning:** The implants are cleaned by the producer and therefore there is no need to clean them or disinfect. If the disinfection is still needed on unused implant, it is necessary to use such disinfectant that does not contain any chloride ions and is suitable for steel implants. The implant must not be damaged (scratched) during disinfection.

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