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SUTURE ANCHOR

The suture anchor is a metal implant which serves for anchoring of damaged tissue with the help of a sewing thread to the bone, e.g. by sewing the damaged parts of the articular capsules. Before the implantation it is necessary to supply the suture anchor with a sewing thread. Then a hole is pre-drilled in the bone and the implant is inserted. A simple set of instruments is used to insert the implant. Various numbers of suture anchors should be used in accordance with the size of the tissue that is being attached.

This implant is made of a titanium alloy Ti 6Al 4V (a material complying with the ISO 5832-3 norm). The product is finished by Type III anodizing, usually of golden color.

Indications

anchoring of soft tissue to the bone:

- refixation of the articular capsule of the shoulder joint, labrum
- hip joint labrum surgery
- others

Contraindications

The surgeon must be well acquainted with the surgical technique.

Conditions that increase the risk of failure:

- an active or anticipated latent infection
- an inflammation in, or close to, the surgery area
- a sensitivity of the patient to the material of the implant (documented or supposed)
- osteoporosis, which does not allow solid anchoring
- obesity or excessive body weight of the patient
- mental or neuromuscular disorders, which can create unforeseen risks of the fixation, especially in the post-operative period
- decreased vascularity, which could impair the blood circulation in the surgery area

Surgical technique

- 1) Find a proper place for the insertion of the anchor
 - nerves, blood vessels, muscle insertions and other similar tissues must not be damaged
 - it is necessary that the bone in the place of insertion is of good quality
 - the place of insertion should be chosen so that the anchor is placed at least partially in the cortical bone
- 2) Perform the incision of the soft tissues
- 3) Use the enclosed drill (diameter 4 mm) to pre-drill the bone into the required depth
- 4) Put (screw) the anchor supplied with a sewing thread onto the Introducer and insert into the pre-drilled hole into the required depth
- 5) Fasten the sewing thread to the part that is being attached, e.g. the labrum of the articular capsule
- 6) Move the anchor deeper into the pre-drilled hole using the Introducing Rod; this will tighten the sewing thread and anchor the tissue to the bone more firmly.

The suture anchor works on the basis of its own self-hooking effect. This means that after inserting it and by applying force to the sewing thread the self-hooking effect takes place and the anchor's movement along the axis in the hole (extraction of the anchor) is prevented. The extent of the resistance against the extraction depends mainly on the quality of the bone. The self-hooking effect is achieved by threading the sewing thread into an eccentrically located opening at the anchor's body. When a force is applied to the thread, the anchor will turn in the hole (anchor gets stuck), and thus the anchor will hook. The anchor is designed in a way that enables the self-hooking effect and maximizes the resistance against an extraction.

Re-implantation of the suture anchor is not expected. In spite of this it is possible to remove the suture anchor with minor difficulties. The Introducer is inserted into the cleaned opening, it is screwed on the thread of the anchor (the Introducer must be tilted in accordance with the turn of the anchor); the anchor is inserted slightly into the hole and aligned, and then gently pulled out of the hole without applying excessive force.

Recommendation: **sewing thread size 2**

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