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Sledujte datum EXPIRACE.

BARBED BROACH – STERILE**WARNING**

1. Products are delivered onto the market in sterile packaging.
2. Before using instruments, check the label for the indication STERILE, the date of EXPIRATION and visually check the integrity of the package. If the package is damaged, return the package to the producer to be repacked.
3. Attention, endodontic instruments are sharp.
4. The root canal instrument may cause breaking of the package during liquidation – when sorting wastes, they must be placed into firm packages – plastic containers or containers secured against opening.
5. Only dentists may work with the instruments.
6. The producer is not liable for any damage caused to the health of the patient or property.
7. Instruments contain nickel, which may cause an allergic reaction in patients sensitive to this metal.

Restriction of repeated processing: Instruments Barbed Broach are designated for single use.

The producer does not recommend repeated use of the instrument due to the possibility of contamination or the transfer of infection. Consequent cleaning of the instrument is not sufficient to remove any possible contamination of the product.

INSTRUCTIONS**Storage and transport:**

Store products in individual packages (blisters) in places where there is no risk of damage to the covering foil or paper. The best storage method is in group packages in paper folders.

Handle the products so as to prevent damage to the package to ensure the sterility of the product. If during handling the package is damaged, return the package to the producer for re-packing.

Monitor the EXPIRATION date.

Checking and testing the function:

Before use, carefully visually check the product and separate all damaged, bent, broken or imperfect instruments.

Storage:

The stored products must be visibly and properly marked to prevent their improper use.

The producer recommends storing the products in paper boxes, in group packages.

Products in sterile packages must be stored so as to prevent package breakages, in particular damage to the covering foil or paper.

When issuing from the store, products must be taken first of all which have been stored for the longest time.

Monitor the EXPIRATION dates.

During the storage of products, it is forbidden:

- to expose products to direct solar radiation, influence of the atmosphere, thermal radiation and mechanical damage
- to use a storage place with a temperature less than 10 °C or more than 40 °C
- to store products in an open area

Other information:

The company MEDIN, a.s. has validated the above-mentioned procedures of cleaning, disinfection, drying and sterilizing. By this validation, these procedures were acknowledged to be satisfactory.

The user is liable for ensuring the necessary materials, and the equipping and training of employees, in accordance with these regulations to achieve the required result.

This *PL0198 Barbed Broach – sterile* is the amendment of *PL0039 Dental rotary instruments – general information*.

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