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Depth gauge

PRODUCT DESCRIPTION


Angularly stable plate may have one of the three basic types of holes for plate fixing to the bone with the help of screws, namely in their optional combination. However, it always contains at least one locking hole.

Locking hole

Self-tapping or self-drilling screws can be introduced into the bone and „locked“ by fastening of the conical thread on the screw head in the same screw in the plate hole.


Such connection is towards the plate angularly stable, the whole system works on the principle of inner fixator. This solution helps to prevent from following problems:

- primary loss of reposition
 - secondary loss of reposition, in the first place in cases of fragmentation fractures without sufficient bone support, low quality or osteoporotic bone
 - periosteum compression and following worsening of blood supply of cortical bone
- Use locking screws up to 3.5 mm in this hole. It is also possible to use standard screw with ball head with stem up to 3.5 mm.

Standard oval hole

Standard oval hole enables to the surgeon to use compression screws for fastening of broken fragments and fracture fixation by axial compression. The screws can be introduced under different angles, lengthwise or transversally. Do not insert locking screws in this hole!

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DISTAL MEDIAL TIBIAL PLATES

Caution: Only a specialist can work with the implants!

Plates

Distal medial tibial plates

Locking screws

Locking screw Ø 2,7 mm self-tapping

Locking screw Ø 3,5 mm self-tapping

Other implants

Locking spacer 3,5 mm

Sleeves

Guide sleeve locking Ø 2,9 mm, L=60 mm

Guide sleeve for rods Ø 2 mm, L=60 mm

Guide sleeve for rods Ø 1,5 mm, L=60 mm


Drills and wires

Drill for small fragments Ø 2,9 mm, L=190 mm

MEDIN K-wire with shank Ø 2 mm, L=300 mm

Autocompression locking hole


Autocompression hole enables dynamic compression in lengthwise axis of the plate at eccentric screw placing. Holes are not regular, compression is always effective in one direction only, from that reason the holes are placed in a way so that the compression would work towards the plate centre.

When using locking screws the hole behaves in the same way as locking hole.

System of angularly stable plates

The system is based on the principle of the inner fixator. By means of the tight connection „screw - plate“, the angularly stable connection is achieved between the plate system and the bone. This transfers part of the weight from the screw-bone connection to the screw-plate connection. Thanks to that it is possible to use the system successfully even in a less quality and an osteoporotic bone.

If the angularly stable plate is used with locking screws, the plate and the screws make tight system together. The plate does not have to be pushed against the bone therefore the bone still gets sufficient supply of blood. Locking screws can be introduced monocortically without enfeebling the strength of the configuration or the fracture fixation, because the screws are firmly anchored in the plate.

PRINCIPLES OF FIXATION

1. Only compression screws

It is about a commonly used method of fracture fixation with the help of plate and compression screws. Fragments of a bone are fastened to the plate, the plate is closely pushed against the bone and stability of the whole system is reached by bone pressure on the bottom part of the plate with the help of compression screws. Shape of the plate is crucial for the correct function of osteosynthesis.

2. Only locking screws

If you use only locking screws then you make a system that works on a very similar principle as outer fixator. When using only locking screws, the plate shape is not what matters but high-quality fracture reposition before tightening of the locking screws.

3. Combination of locking and compression screws

- Use the plate on the base of the common plating principle; first perform fracture reposition and fixation by compression screws. (You can use the spacer to keep the space between the bone and the plate). Then screw the locking screws to reach the angularly stable plate fixation towards the bone.
- After the fracture reposition and plate fixation by locking screws, it is possible to fasten the separated fragment (e.g. multi-fragmentation or splintery fractures) to the fixed part of the bone by compression screw.
- After fixation of one fracture fragment by locking screws, it is possible to push the second fragment to the first one with compression screws and eventually after that to fix even the second fragment with other locking screws.

Note: Do not use locking screws closely to both sides of the fracture because the strain, to which the plate is exposed, would be concentrated only to the space between these screws and the plate could get broken. It is generally recommended to leave from 1 to 3 holes empty above the fracture area.

Note: To use the locking plate system correctly, it is necessary to understand the principles of fixation and their biomechanical properties. Draw your attention to pre-operation planning.

INDICATION

The plate is used for fixation of distal tibia fracture parts. Angularly stable screws are mainly used at articular fractures. The angles of the screws introductions are chosen in a way so that they would fix as many fragments as possible and they would minimise the possibility of screws pulling out of the bone. The distal part of the plate is the same at all the sizes, there are 10 holes for locking screws. The straight shank part has from 2 to 10 holes and 1 oval hole. The longer plates enable fixation in the area of diaphysis.

CONTRAINDICATION

See general information about implants using – Instructions nr. PL0088.

Combination with other products

It is possible to combine plates with screws and wires of MEDIN, a.s., providing the instructions and producer's recommendations are followed, particularly compatibility of the materials. (it is not allowed to combine products of different materials). The products are not determined to be combined with products of other producers.

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
PLATE INTRODUCTION AND LOCKING WITH LOCKING SCREWS

Despite of the fact that locking screws are introduced bicortically in most cases, it is possible to introduce them even monocortically in quality and strong bones. Monocortical introduction of the screws is mainly used in case when the screw is not desired to go through the whole bone e.g. in the area of epiphysis.

Pay attention to locking screws, they do not behave as standard compression screws. It is necessary to realise that during fastening the screw and locking the screw's head in the plate hole, the screw is strongly fixed and the force compression of bone fragments does not happen. Therefore, careful anatomical bone reconstruction is necessary (mainly in intra-articular fractures) with the help of compression (cortical or spongy) screws, K-wires or cannulated screws. If the intra-fragment compression is required, using of standard screws is also necessary.

1. Reposition the fracture

Reposition the fracture under X-ray control. If essential, fixate it with K-wires MEDIN or forceps.

2. Insert the spacer

To minimise the contact of the plate with the bone, it is possible to place the spacer in the hole with thread before fixation of standard screws, it ensures the space between the plate and the bone 2 mm. It is possible to remove it again after plate securing with locking screws.

3. Insert guide sleeve locking

Screw the guide sleeve locking carefully into the thread of chosen plate hole until it is tightly hold. Use the guide sleeves locking for manipulation and correct plate placing.

4. Place the plate and fix it temporarily

Place the plate where it should be with the help of guide sleeves locking and fix it with K-wires MEDIN or forceps.

Correct position of the screws is necessary for right function, mainly screw alignment with the plate hole. To determine correct screw position you can use K-wires. Screw locking sleeve 2.9 mm into required plate hole and insert guide sleeve for rods 1.5 mm into it and introduce the wire. Use RTG to picture the wires – they represent final position of locking screws.

5. Pre-drill the hole

Remove the guide sleeve for rods Ø 1.5 mm and leave the drilling sleeve Ø 2.9 mm in the place. With the help of this sleeve, pre-drill the hole for screw. Drill only until first cortex then you can screw the screw easily thanks to the self-tapping thread on the tip.

Carefully with a drill (Ø 2.9 mm for screw of Ø 3.5 mm; for screws Ø 2.7 mm pre-drill the bone

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with wire of Ø 2 mm and using guide sleeve) pre-drill hole in the bone with the help of attached guide sleeve. Use the electric drill for mechanical pre-drilling. Use the hand-handle of the drill for more exact control above depth of drilling and to prevent the soft tissues from damage on the opposite side of the bone.


6. Measure the hole depth and choose the correct screws

Unscrew the guide sleeve carefully. In case that it is not possible to unscrew the sleeve only by hand from some reason, you can help yourself with guide sleeve for wires – insert it into the hole of the drilling sleeve handle and use it as an elevator. Use the wire hole for screw depth determination with gauge.

7. Insert the screw

Before introduction of the first locking screw, make sure that the plate is firmly fixed - when fastening the screw the rotation could happen – the plate round the screw axis. Make the screws tight with the help of hexagonal moment screwdriver.

Pre-drill the holes for locking screws with the help of the guide sleeve. For proper plate function, it is necessary to introduce the screws in the plate in a way that the lengthwise axis of the screw is identical to the axis of the corresponding hole in the plate. At a divergence smaller than 5° the connection of the screw's head and the plate is decreased to 70 %.

8. Make sure everything is all right

Check that the locking screws are tight. Remove instruments necessary for temporary plate fixation.

9. Plate removal

When removing the plate, loose all the screws first, after that start to take them out one after the other. In this way, you prevent from possible rotation of the plate round the axis of the screw that was taken out as the last one.

If the hexagon head of the screw gets stripped or during other problems with screws removing, use the instrumentation for removal of MEDIN screws.

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