## HŘEB CALCANEUS C-NAIL

### ÚČEL POUŽITÍ

#### Určený účel použití

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### PURPOSE

#### Intended use

C-NAIL is an intramedullary nail used for fixation of a fractured heel bone, or calcaneus. The principle is to stabilize a fragment of a fractured calcaneus in conjunction with locking cross screws, which creates angular stable fixation. Introducing the nail enables a mini-invasive procedure. Fixation is possible with up to seven screws. Using this operating technique can resolve both simple and complicated fragmentation of intraarticular fractures.

#### Description

C-NAIL is produced from a titanium alloy or implant steel. The design of the nail is separated of the left leg (indicated by L) and the right leg (indicated by R). The difference between individual designs is found in the differing direction of holes for screws into the sustentaculum. The nail length is 65 mm with a diameter of 8 mm. It is fitted with seven holes for screws:

- 2× sustentacular screws
- 2× superior screws
- 3× lateral screws

The nail is used in a conjunction with special screws with a flat head with a diameter of 4 mm or 3.5 mm with low-profile thread. The operator chooses the length of the screws based on location. The proximal end of the nail is closed with a stopper after introducing the nail and screws.

Special instrumentation is applied to introduce the nail and screw. FOLLOW INSTRUCTIONS FOR USE OF THE INSTRUMENTATION AS LISTED IN THE PACKAGE LEAFLET. Prior to use, assure proper configuration of the aiming device arms based on type of applied nail. Prior to use, check for correct aiming in the configuration of the aiming device, nail and guide sleeves. The drill running through the drill sleeves must move freely through the holes in the nail. In the event of damage or any inconsistencies, it is not possible to use the implant for surgery.

#### Principle of fixation

After reduction, individual fragments of the calcaneus are fixed with a nail and locking screws.

Detailed use of the implant and instrumentation during surgery is described in the operating instructions.

#### Indication

Extraarticular and intraarticular fractures of the calcaneus according to Sanders I–IV, Joint Depression type, Tongue type.

#### Contraindication

1. Insufficient quantity or quality of bone, which may prevent proper fixation of the bone.
2. Any fully developed or putative latent infection.
3. Previous course of infections.
4. Any mental illness or neuromuscular illness of the patient that would lead in the postoperative phase to unacceptable risk in relation to unsuccessful fixation or complications.
5. Reduced vascularity, which would prevent the necessary blood supply to the fracture or operated area.
6. Malignant growth at the point of fracture.
7. All cases where the implant would come in conflict with anatomical structures or physiological functions.
8. Use of the steel implant if the patient is allergic to Ni.
9. The calcaneus shorter than 65 mm.

#### Absolute:

1. Skin or bone infection.
2. Malignancies.
3. For children with unclosed calcaneal sulcus.

#### Relative:

1. Diabetes mellitus.
2. Ischemia of the lower extremities and other circulation disorders of the lower extremities.
3. Uncooperative patient.

### Surgical procedure

#### 1. Preoperative planning and principle of surgical technique

The underlying principle of the operation is fixation of a fractured calcaneus with a nail and locking or angular stable screws. From a small subtalar access point after visualization of the fractured posterior joint surface, perform reduction and fixation using individual screws while under X-ray control. After restoring the Böhler angle, fix the fragments temporarily with K-wires. After introducing the aiming wire to the marrow cavity of the calcaneus under the attachment of the Achilles tendon towards the center of the calcaneocuboid joint, you then drill a hole for the nail. Fix the introduced C-Nail under X-ray control with the help of X-RAY transparent arms first to the sustentacular fragment, consequently fix this with both lateral and superior screws.

Within the framework of preoperative planning, you must verify the appropriateness of this surgical technique. One of the parameters that must be verified is the length of the calcaneus, which may not be shorter than 65 mm. In the opposite case, it is not appropriate to introduce the nail due to non-standard length. We recommend taking X-rays in the lateral, axial, dorsoplantar and Broden projections, and CT of the calcaneus. It is appropriate to perform a comparative lateral X-ray of the healthy extremity. It is not appropriate to apply on patients with unclosed calcaneal sulcus.

#### 2. Preparation

##### 2.1 Preparing instruments

###### 2.1.1 Preparation of aiming device

On the basic body of the aiming device, aiming device arms are set into a position based on the chosen side of the operated extremity. The aiming device arms are fixed with fixation screws, which are tightened reasonably to the basic body with the help of tightening rods.

##### 2.2 Assembling of the aiming device and introducing the nail

###### 2.2.1 Assembling the aiming device with nail

The nail is placed in the aiming device grooves, and is firmly fixed using the aiming device screw.

###### 2.2.2 Assembling the aiming device with the sleeves and inspecting the aiming device with nail

Prior to use, inspection must be performed of the position of the aiming device arms against the tight nail. Inspection is performed for all seven holes.

##### 2.3 Preparing the implant

###### 2.3.1 Choice of nail

Nails are constructed separately for the left and right leg, and are not supplied sterile. Before use, nails must be disinfected, washed and sterilized.

###### 2.3.2 Preparation of instruments

Prior to surgery, it is necessary to check the completeness and function of all instruments, which are supplied in unsterilized condition. Before use, nails must be disinfected, washed and sterilized.

#### 3. Operating technique

##### 3.1 Position of patient

The patient lies on his or her side so that the operated extremity would be slightly bent in the upper position. The operated extremity is supported by a soft cylindrical pillow from the medial side. The other extremity is bent at the knee to not interfere with using X-ray equipment.

##### 3.2 Drawing the position of the bone and tendons

The drawing is performed from the lateral side for performing access to the subtalar joint.

##### 3.3 Incision for reduction of the posterior joint surface

For reduction of around 3–4 cm long incision under the attachment of the external ankle to the base of the fifth metatarsal bone. When making the incision, exercise extreme caution to keep from damaging the peroneal tendons.

##### 3.4 Reduction of the posterior joint surface

Reduction of the rear joint surface using an elevating device, fixation of reduced fragments. Configuration of the joint space is recommended from the medial fracture towards the lateral fractures.

##### 3.4.1 Additional fixation using K-wires

The current fixation of reduced fragments of the rear joint space is secured using K-wires. Correctness of the congruence of the joint space is now checked in the Broden projection.

##### 3.4.2 Fixation using cancellous screws

Fixation of the rear joint space by one or two separate cancellous screws with diameter of 4 mm introduced under the joint surface.

##### 3.5 Leveling the Böhler angle and axial deviation (varosity or valgosity)

###### 3.5.1 Inserting the Schanz screw to the tuber calcanei

Using the torque of the inserted Schanz screw, you can correct the varose or valgose position or shortening of the calcaneus, and restore the Böhler angle.

###### 3.5.2 Fixation of the K-wire position

K-wires must not prevent consequent introduction of the nail. The recommended fixation is plantar into the talus using a K-wire. If necessary, use of another K-wire is possible, introduced in parallel direction with the calcaneus.

##### 3.6 Introducing the nail

###### 3.6.1 Incision

The incision at around 15 mm in length is performed vertically under the attachment of the Achilles tendon.

###### 3.6.2 Drilling the guide wire for the nail hole

The direction of drilling the guide wire is determined by the center of the calcaneocuboid joint. Drilling is performed using a guide wire through a protective sleeve, in which an insert is placed with a hole for the guide wire. Inspection of drilling and drill depth is performed in two projections by X-ray amplifier. The guide wire is inserted to the subchondral area of the calcaneocuboid joint.

###### 3.6.3 Drilling the hole for the nail

Drilling is performed through a hollow guide wire cutter 8 mm in diameter. The soft tissues are protected during drilling by the guide sleeve. The depth of drilling is measured at the X-ray side projection approximately 5 mm in front of the CC joint. It is possible to calculate the drilling depth on the cutter scale against the guide sleeve. After drilling the hole, the cutter and guide wire are extracted.

###### 3.6.4 Introducing the nail with the aiming device

The nail with the aiming device is introduced to the drilled hole. This is performed using light axial pressure by rocking axial motion of the aiming device with the nail.

###### 3.6.5 Adjusting the nail position and introducing K-wires to the sustentaculum

The correct nail position is given by the direction of the aiming device wires to the sustentacular fragment. Prior to introducing the temporary wires to the sustentacular fragment, it is necessary to correctly set the depth and rotation position of the aiming device with the nail. First the temporary wire is drilled through the proximal hole on the sustentacular arm through the incision enabling introduction of the guide sleeve (When drilling all holes, it is necessary to make sure the drill sleeve is introduced to the incision together with a conical trocar for parting the soft tissues. Prior to inserting the K-wire or drill, the protective guide sleeve must be resting firmly against the bone). After drilling the hole, check the position and depth of the temporary K-wire by X-ray amplifier in two levels (lateral and axial). If the position is not correct, it is necessary to remove the wire and perform correction of the aiming device position. Next, the second hold is drilled for the distal screw into the sustentaculum tali. Depth of drilling is inspected using the X-ray amplifier. The fixation screw is introduced into the hole. The screw length is measured using a calibration drill or gauge by subtracting the length from the edge of the guide sleeve. Introducing the screw is checked in two projections under X-ray amplifier.

###### 3.6.6 Introducing the remaining screws

The other fixation screws are then introduced. During introduction, it is necessary to carefully measure the length of necessary screws so that they are introduced bicortically. Screws should not protrude through the second cortex in order to minimize irritation of soft tissue. The screws are specially designed with a flat head to prevent their undesired embedding into the bone. Screws are placed on the screwdriver and introduced to the drilled hole through the guide sleeve (blue strip). The screw head rests upon the cortex by gently tightening with the screwdriver. In the correct position of the sleeve resting against the cortex, the screw depth is given by the scale on the screwdriver handle. A hand-held screwdriver is used exclusively for introducing screws. Choice of using all holes is up to the surgeon based on the indication. It is recommend to use all screws for solid angle fixation of the fragments. When introducing the center screw to the lateral side, it is necessary to exercise caution in order to prevent any damage to the peroneal tendons.

###### 3.6.7 Inspecting fixation

After introducing all screws, the fixation is checked as well as the length and position of screws in the lateral, axial, dorsoplantar and Broden projections. Next, the temporary fixation wires are removed.

###### 3.6.8 Removing the aiming device

The aiming device is released from the nail by loosening the clamping screw. A rod inserted in the screw head hole is used for loosening.

### 3.6.9 Introducing stoppers

A stopper is introduced to the hole in the proximal part of the nail. The length of the stopper is based by the nail as it is embedded (0 mm; 5 mm; 10 mm; 15 mm; 20 mm). The stopper should be flushed with the outer cortex.

### 3.7 Final inspection of reposition, position and length of screws in the lateral, axial and Broden projection.

### 3.8 Completion of the operation

The wounds are rinsed with a sufficient quantity of physiological solution, or a Redon drainage catheter is introduced. Closing the operating access point is performed in individual layers. The final step is application of a sterile dressing, plaster splint or elastic bandage. The chosen variation depends on the decision of the treating physician based on stability of fixation of the fracture and condition of soft tissue.

### Post-operative care

The treating physician determines a rehabilitation plan based on results of the operation and consequent checkups. It is recommended to commence rehabilitation immediately on the first day after operation or after removal of the plaster splint. It is always necessary to respect the fact that the applied osteosynthetic material is only for fixation and is not designed for receiving the full weight of the patient. Full weight is recommended only after creation of a strong bony callus, so that the bone tissue that has grown over receives the full weight. If creation of a sufficiently strong bony callus does not occur, it is not possible to allow full weight on the extremity due to the danger of fixation failure and damage to the osteosynthetic material, and it is necessary to choose an osteosynthesis method based on assessment by the treating physician.

Magnet therapy is recommended to accelerate healing.

### Removing osteosynthetic material

1. It is recommended to extract the implants. All osteosynthetic material however can be left in the patient if there are no reasons for its removal, such as infection, impingement syndrome, pain, etc.
2. Extraction procedure:
  - Loosen stopper.
  - Loosen fixation screws.
  - Nail extraction.

### WARNING

**The aiming device arms are produced from an X-ray-transparent material. Sterilization can be performed only by using steam at a max. temperature of 134 °C, overpressure of 304 kPa for a period of 10 minutes!**

**When applying higher temperatures, permanent deformation and loss of function of the arms may occur!**

When cleaning and disinfecting, it is necessary to avoid agents containing chloride ions, which may cause corrosion of the steel parts.

Implants – different types of materials must not be combined (steel-titanium) – incorrect  
(titanium-titanium), (steel-steel) – correct

Products may not be combined with products of different manufacturers.

### READ CAREFULLY!

### Preventative measures:

1. The implant must never be used twice! Previous stress may decrease its fatigue limit, and if applied twice, it may lead to failure of osteosynthesis.
2. The implant is only for fixation and may not receive the full weight of the patient. After healing and fixed bone fusion, it is recommended to remove the implant. Extraction depends on the judgment of the physician providing treatment.
3. The operator must be familiar with the surgical procedure while using prescribed instrumentation. Only a person who has undergone instruction in accordance with this instruction manual may use these medical devices.
4. Only use the prescribed instrumentation for performing the operation.
5. Consider the appropriateness of using the implant in patients where it is not certain that they will cooperate with the physician after the operation (ex. alcoholics, addicts, the mentally ill).
6. Protect the implant from deformation and scratching, when handling, avoid contact with other metal materials

and chemicals.

7. Do not use implants with a given patient that are different in terms of material, because contact of various metals may accelerate the corrosion process, which may lead prematurely to a decrease in the effectiveness and service life of the implant, as well as to an increased quantity of metallic compounds that are released in the body of the patient and could evoke a toxic reaction.
8. Decide on the possible premature removal of the implant if you suspect poor healing of the bone or that it is becoming loose, bending, or breaking in the bone, because this could lead to serious injury to the patient.
9. Decide on removing the implant for internal fixation while taking into account such factors as the risks for patients in consequence of an additional surgical procedure, and difficulties during and after implant removal.
10. Patients who smoke or use products containing nicotine have increased frequency of the fractures failing to heal. It is necessary to notify such patients of this fact, and warn them of possible consequences.
11. The appearance of implant failure is more probably amongst heavy, physically active persons.
12. The manufacturer bears no liability if the product is used in a manner other than as stated in the instruction manual, or in the surgical procedure.
13. The manufacturer bears no liability if the product is used in combination with implants of other manufacturers.

### Instructions to the patient:

1. Notify the patient that he or she must follow precisely the instructions of the treating physician.
2. Inform the patient on the course of the operation, on potential risks and on the rehabilitation program.
3. Instruct the patient when he or she can apply full body weight on the bone with internal fixation, in order to prevent damage to the implant.
4. Instruct the patient to maintain a reasonable level of activity in order to keep from refracturing the calcaneus.
5. Inform the patient when probable removal of the implant is to occur.
6. Inform the patient to maintain a reasonable level of activity after implant removal, because refracturing could occur before the bone completely heals over and before the bone cavity becomes full.

### Possible side effects:

1. Poor fusion of fragments, suppressed revascularization and poor formation of bony callus may cause loosening, bending or breaking of the implant or premature loss of solid bone fixation.
2. Early or late infection, both deep and superficial.
3. Avascular necrosis.
4. Damage to nerves and tendons as a result of the surgical procedure.
5. Allergy to Ni.
6. Development of arthrosis.
7. A risk of incorrect diagnosis when the diagnostic devices CT or MRI are used where the presence of steel implant in the given locality causes a wrong image.

Serious complications may occur with every surgical operation when removing the results of multiple injuries. These complications include urogenital, gastrointestinal and vascular disorders including thrombosis, bronchiopulmonary disorders including embolism, myocardial infarct and death.

**Limitation of repeated use:** The implants are designed for one-time use in humans. If the implants were prepared for a surgical procedure and became contaminated, they may be prepared for use again through a process of cleaning, disinfection and sterilization described hereunder.

### INSTRUCTIONS

#### NOTICE:

The implant is supplied UNSTERILE. It is wrapped in polyethylene packaging. It is cleaned and disinfected by the manufacturer. Prior to first use, remove the implant from the packaging and perform sterilization based on the regulation. If the packaging is not fully intact, it is necessary to return the implant to the manufacturer for repackaging.

**Place of use:** Specialized facility for health-care providers, mainly surgical and orthopedic facilities.

**Storage and transport:** Protect the implant from deformation and scratching, when handling, avoid contact with other metal materials and chemicals.

### Preparation for cleaning:

#### IMPLANTS

Implants are cleaned by the manufacturer, so it is not necessary to clean or disinfect them. If it becomes necessary to disinfect an unused implant, it is necessary to use a disinfectant containing no chloride ions and which is appropriate for use on steel instruments. Damage (scratching) of the implant must not occur during disinfection.

#### INSTRUMENTATION

1. Remove the arms from the aiming device body.
2. Perform mechanical cleaning with water and brush.
3. Rinse the instruments with pressurized water.
4. Submerge the instruments in a disinfecting solution for 20 to 30 minutes. The recommended disinfection agent is *Sekusept pulver, NeodisherseptoMED*. Danger of damaging instruments exists if a different solution is applied.
5. Rinse again with pressurized water until only clean water runs off. Rinsing may be performed in a conventionally supplied pressurized washer for connecting instruments to the jets.
6. Blow dry the instruments with pressurized air.

**Disinfection:** We recommend either a washing appliance and thermal or thermochemical disinfection, or hand washing of the implants and instruments after chemical disinfection with an agent having virucidal effects, whereas it is possible to join both stages by use of disinfecting agents with combined washing effect. It is possible to use the disinfecting agent in accordance with the instructions listed on the label of this agent. When preparing the disinfectant and washing solutions, you must uphold the procedure given by the manufacturer. Instruments designed for disassembly are disassembled and each part is then considered to be separate. After performing disinfection, the implants and instruments are rinsed with pure demineralized water to remove chemical residue remaining on the implants and instruments.  
Recommended agents: *Sekusept pulver, NeodisherseptoMED*

**Drying:** If drying is a part of a machine washing and disinfection cycle, the temperature may not exceed 120 °C.

**Maintenance:** After performing each operation, it is necessary to check individual parts of the instrumentation. It is also necessary to check cutting instruments for bluntness. In the event of damage, it is necessary to send the parts to the manufacturer for performing repairs. After 10–20 performed operations, it is recommended to contact the manufacturer for performing a service inspection.

**Inspecting and testing function:** Prior to using the implant, it is necessary to visually inspect for completeness and proper function.

**Sterilization:** If not stated otherwise, sterilization is recommended with moist heat in a steam sterilizer equipped with an antibacterial filter at a temperature of 121 °C, overpressure of 205 kPa for a period of 20 minutes or a temperature of 134 °C, overpressure 304 kPa and for a period of 10 minutes.

**Further information:** The company MEDIN, a.s., validated the above-mentioned procedures of cleaning, disinfection, drying and sterilization. By this validation, these procedures have been recognized as qualified for preparation of medical devices designed for repeated use.

The user is responsible for acquiring the necessary materials, equipment and training of employees in accordance with these prescribed procedures, in order to achieve the required result.

The user may also apply other procedures, but the effectiveness of these alternative procedures must be verified by a qualified laboratory.

**Storage:** Store the packed products in a dry and clean environment without extreme temperatures and humidity, out of the reach of direct sunlight.

**Liquidation of the implant:** The implant is considered hazardous waste after use. The user is responsible for performing measures assuring safe handling and liquidation of the product.  
After drying, used implants are discarded and put aside as potentially hazardous waste.

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