## VAROVÁNÍ

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## CANNULATED SELF-TAPPING SCREWS HA 3.5xLxK mm

## INTENDED USE

## Description of the implant

Cannulated screws HA 3.5 are not fully threaded bone screws intended for use in bone surgery and traumatology. They are made of implant steel (under ISO 5832-1) or Ti6Al4V alloy (under ISO 5832-3). The screws are provided with a channel that allows guidance of these screws with K-wires. Opening for the guide wire corresponds to the Ø 1.2 mm K-wire, resp. Ø 1 mm K-wire. These types of screws are not fully threaded, the thread length increases with the length of the screw (see the catalog of MEDIN, a.s.). Raised head is equipped with 2.5 mm socket.

## Compatibility

Cannulated screws HA 3.5 are compatible with K-wire of 1.2 mm and K-wire of 1 mm. These screws are intended for connection with bone plates or separately.

## PL for K-wires:

PL 0109 – Kirschner wires, K-wires

## Indication

Screws are intended for simple fractures of the upper limb, ankle and foot.

## Contraindication

- Insufficient quantity or quality of bone, which may prevent correct fixation of bone.
- Any fully developed or presumed latent infection.
- Previous course of infections.
- Any mental disorder or neuromuscular disease of the patient, which may result in unacceptable risk related to fixation failure or complications in the postoperative period.
- Impaired vascularity, which may prevent the necessary blood supply to the fracture or surgery site.
- Malignant growths in the fracture site.
- All cases where the product is in conflict with anatomical structures or physiological functions.
- Using a steel implant if the patient is allergic to nickel.

## Absolute:

- Skin or bone infection.
- Malignancies.

## Relative:

- Diabetes mellitus.
- Uncooperative patient.

In cases of implantation in the lower limbs, a relatively low bearing capacity of the screw due to strain applied to the fracture and prescription of appropriate rehabilitation program must be taken into account.

## Fixing principles

Perform reduction of fragments into the anatomical position and fix it, according to the size of the fragment, e.g. with two K-wires. Insert a guiding Ø 1.2 mm K-wire. Check the correct position of the wire with X-ray, K-wire position corresponds to the final position of the screw. Measure the length of the screw with another equally long K-wire or depth indicator. Pre-drill the hole with a Ø 2.4 mm drill and insert the screw using a cannulated screwdriver. After the insertion and tightening of the screw you can compress with further tightening. But only if the threaded part of the screw is fully inserted into the distant fragment and the non-threaded part is introduced into the closer fragment. To prevent rotation of the fragment, it must be fixed with a second screw, resp. another implant, which will function as anti-rotation.

Screws may also be used in combination with plates, or as transfixation material, before the introduction of the „main“ implant.

## Post-operative care

Depending on the fracture nature, fix the leg with an orthosis or plate in the postoperative period. Recommend optimal rehabilitation program to the patient, based on individual assessment.

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## WARNING

When cleaning and disinfecting, it is necessary to avoid agents containing chloride ions, which can cause corrosion of steel parts.

Implants (plate and screws) from different materials may not be combined

steel x titanium ❌

titanium x titanium or steel x steel ✅

Products may not be combined with products by other manufacturers.

## Reprocessing limitation

Implants are disposables for humane use. If the implants were prepared for surgery and not contaminated, they can be re-prepared for use by cleaning, disinfection and sterilization process as described below.

## INSTRUCTIONS

**WARNING:** The implant is supplied NON-STERILE. It is packaged in a polyethylene bag. It comes cleaned and disinfected from the manufacturer. Before first use, remove the implant from the package and perform sterilization as prescribed. If the package is broken, it is necessary to return the implant to the manufacturer for repackaging.

**Place of use:** Professional workplace of health care providers, namely surgical, trauma and orthopedic departments.

**Storage and transport:** Prevent deformation and scratches to the implant; avoid contact with other metallic materials and chemicals when handling the implant.

**Preparation for cleaning:** Implants are cleaned by the manufacturer, and therefore need not be cleaned or disinfected. Should it be necessary to re-sanitize an unused implant, use a disinfectant that contains no chloride ions and is suitable for steel instruments. The implant must not be damaged (scratched) during disinfection.

**Cleaning – automated:** No special requirements.

**Cleaning – manual:** No special requirements.

**Disinfection:** Either device for washing and thermal or thermochemical disinfection is recommended, or hand washing of implants after chemical disinfection with virucidal effect, while the two stages can be combined using disinfectants with combined washing effect. Disinfectant can be used in accordance with the instructions on the label of the product. In the preparation of disinfectants and cleaning solutions, the procedures specified by the manufacturer must be followed. After disinfection, the implants are rinsed with demineralized water to remove chemical residues remaining on the implant.

*Recommended products:* Sekusept pulver, NeodisherseptMED

**Drying:** No special requirements.

**Maintenance:** No special requirements

**Inspection and testing of function:** Before use, the implant should be inspected visually.

**Packing:** No special requirements.

**Sterilization:** Unless otherwise stated, it is advisable to use steam sterilization in a steam sterilizer equipped with antibacterial filter at 121 °C, pressure 205 kPa for 20 minutes or at 134 °C, pressure 304 kPa for 10 minutes.

**Storage:** The implant must be stored in a dry place without direct exposure to chemicals.

**Implant disposal:** The implant is considered hazardous waste after use. User is responsible for implementing measures of safe handling and disposal of the product.

After drying, used implants are discarded and disposed of as potentially hazardous waste.

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