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PELVIC PLATES STRAIGHT, CURVED AND OMEGA

INTENDED USE

Description

Straight, curved and OMEGA pelvic plates are used to stabilize fractures of pelvic bones. The OMEGA plate is mainly used for reconstruction of the acetabulum quadrilateral area using the Stoppa's approach, or for treatment of fractures of the acetabular anterior column and to stabilize fractures of the pelvic ring. These implants are only available in steel (steel alloy e.g. 1.4441 according to ISO 5832-1)

Compatibility

The plates are used with 3.5 mm HA cortical screws or 4 mm HB cancellous screws. Implantation of the plate requires special instrumentation and instrumentation for insertion of the above screws. Pelvic plates must be fixed with screws; pelvic plates can also be mutually combined. Use PLO139 to insert screws.

Operation procedure

Straight and curved plates

All pelvic joints (acetabular, sacroiliac) to be connected must be properly reduced. Plates must be properly selected according to the desired curvature, or re-shaped with special bending pliers to fit snugly the area of operated bones; at the same time, alternating bending is to be avoided anywhere on the plate, as they reduce the ultimate strength of the material (cyclic fatigue). Plates can be re-shaped in the thickness plane (ripple following the surface of the bone) and in the width plane, if basic semi-diameters are not sufficient (flat and curved plate with a midline semi-diameter of 88 mm) and spatially (in the shape of the screw area). When re-shaping in the thickness plane, deformation at the screw holes and among them, i.e. at the notches, must be avoided. When fastening the plate with screws, moving fragments to the plate must be avoided to prevent their dislocation. The type of screw depends on the type of bone to be fixed (HA 3.5 cortical screw or HB 4 cancellous screw), screw length depends on the thickness of the bone in the direction of the screw (must not be too large because of pelvic viscera).

OMEQA plates

The upper arm of pubic bone and the quadrilateral area of acetabulum are released from Stoppa's approach. The plate is placed under the rounded linea arcuata. Extended part of the plate is supported by the quadrilateral area and stabilizes the central dislocation of the hip.

The plate must be shaped properly to fit the shape of the bones at the site of attachment, to follow the shape of bone mass as much as possible. When shaping the plate, it is advisable to use the screw holes. Bends of the plate should be continuous, without breaks (e.g. at the holes).

When shaping the plate, it is inappropriate (forbidden) to bend it back. There is a danger of microcracks developing in the place of bending, which can subsequently lead to plate breakage due to fatigue during loading.

It is recommended to fix the re-shaped plate with HA 3.5 cortical screws inserted in bicortical way. The holes in the plate allowed the insertion of screws in different directions depending on the desired fixation. Alternatively, it is possible to insert HB 4 mm cancellous screws.

Note:

In order to facilitate bending of the pelvic plate, it is possible to use a tape that can be easily modeled according to the needs by hand and the plate can then be shaped according to this model. The tape, however, must never be used as an implant.

The necessary instrumentation

During the operation, apply common classic surgical instruments to open the operating area, as well as retractors and surgical instruments to reduce bleeding.

Within the implantation of pelvic plate use hooks for pelvic surgery, tape, pliers, bender (pin, fissure and fissure perpendicular) in several designs for plate shaping, bending pliers, ball retainer jabber, 2 mm k-wires, protective sleeve, drill of 2 and 2.7 mm diameter, tap, depth gauge, and 6HR 3.5 mm screwdriver.

Conclusion

Recommended surgical approach, especially for OMEGA plate's surgery, follows Stoppa's approach, which allows a good overview of the operational space (namely fair visibility in the surgical area) with the possibility of quality stabilization of bone fragments. The plate can be fixed with screws close to the fracture, which is the pre-requisite for high-quality and solid osteosynthesis.

Note:

The company MEDIN, a.s., Nové Město na Moravě, offers re-shaping and shortening of OMEGA plate at the hip joint centers according to received CT scan. To use this service, please contact our sales representative for exact conditions.

Indication

- Transverse acetabular fractures
- T-shaped fractures of the acetabulum
- Fractures of the anterior wall and anterior column of the acetabulum
- Complex fractures of the acetabulum
- Fractures of the pelvic ring
- Central dislocation with femoral head medialization
- Fractures of the posterior wall and posterior column of acetabular cup
- Pubis symphysis fractures, rear vertical fractures of the hip and pelvic fractures in the area adjacent to the sacrum

Contraindication

- Insufficient quantity or quality of bone, which may prevent correct fixation of bone.
- Any fully developed or purported latent infection in the surgical area or in its vicinity.
- Previous course of infections.
- Any mental disorder or neuromuscular disease of the patient, which may result in unacceptable risk related to fixation failure or complications in the postoperative period.
- Reduced bone stock, impaired by disease, infection or prior implantation, which cannot provide proper support for implants
- Impaired vascularity, which may prevent the necessary blood supply to the fracture or surgery site.
- Malignant growths in the fracture site.
- All cases where the medical product is in conflict with anatomical structures or physiological functions.
- Using a steel implant if the patient is allergic to nickel (documented or suspected allergy).
- Obesity or overweight patient

Post-operative care

Recommend optimal rehabilitation program to the patient, based on individual assessment.

WARNING

	Implants (plate and screws) from different materials may not be combined
steel x titanium 	titanium x titanium or steel x steel 
Products may not be combined with products by other manufacturers.	

Preventive measures:

- The implant should never be re-used! Previous stress can reduce the fatigue strength and the re-implantation may cause failure of osteosynthesis. Used implant may be contaminated and its re-use can cause an inflammatory response of the body. It can also serve as a media of disease of the first patient.
- The implant cannot indefinitely withstand the level of activity corresponding to a normal healthy bone. It can be used until the bone is healed and afterwards, in active younger patients the implant for fixation of fractures is usually removed. In older patients the surgeon can consider whether to retain the implant for fracture fixation in place after healing or risk additional surgery for its removal, even if the implant has lost its effectiveness.
- The surgeon must be familiar with the operation procedure for the specified instrumentation.
- Use only the specified instrumentation for the surgery; to work well, it must not be damaged.
- Consider the use of implants in patients where there is no certainty that they will cooperate with their physician after the surgery (such as alcoholics, drug addicts, and the mentally ill).
- Prevent deformation and scratches to the implant; avoid contact with other metallic materials and chemicals when handling the implant.
- If it is necessary to change the shape of the implant, it can only be performed as specified.
- Do not use implants from different materials in one patient because contact between different metals may accelerate the corrosion process, which may lead to a premature reduction in the efficiency and life of the implant, as well as an increased amount of metal compounds released in the body of the patient, which may cause toxic reactions.
- Opt for early implant removal if poor bone healing, implant loosening, bending or breaking in the bone is suspected, as it can result in serious injury to the patient.

- Make sure to keep all shapes and sizes of implants that are likely to be used them available at the medical facility.
- Consider additional preventive means which reduce the weight of the load in order to minimize premature loading device for fracture fixation, until the X-ray visible connection bone.
- Opt for implant removal for internal fixation, taking into account factors such as the risks to patients due to additional surgery and difficulties during and after implant removal.
- Patients who smoke or use products containing nicotine on regular basis have an increased frequency of unhealed fractures. Such patients should be advised of the fact and warned about the possible consequences. The incidence of implant failure is more likely in heavy, physically active people.
- The manufacturer is not liable if the product is not used in accordance with the instruction manual or operating procedure.
- When using visual methods such as X-ray, CT scan, magnetic and other investigative methods to be followed and follow the instructions of the manufacturer of these devices.
- The manufacturer is not liable if the product is used in combination with implants from other manufacturers.

Instructions for the patient:

- Tell the patient that they must strictly follow instructions of the attending physician.
- Instruct the patient about the course of surgery, risks and rehabilitation program.
- To prevent damage to the implant, instruct the patient when the bone with internal fixation can be fully strained.
- Instruct the patient that the bone is prone to damage during the healing, particularly in places where the implant ends.
- Tell the patient when the implant will be probably removed.
- Instruct the patient to keep to reasonable level of activity after implant removal, as otherwise a re-fracture may occur before the bone completely heals and the bone cavity is filled.

Possible side effects:

- Failure to join the ends of fractures, osteoporosis, suppressed revascularization and poor bone formation may cause loosening, bending or breakage of the implant, as well as premature loss of solid bone fixation.
- Bad connection of the point of fracture resulting from improper alignment.
- Increased fibrous tissue reaction in the fracture site due to unstable shattered fractures.
- Early or late infection, both deep and superficial.
- Avascular necrosis.
- Shortening of the affected bone (fracture) side.
- Subclinical nerve damage may occur as a result of surgical trauma.
- Allergies to Ni
- The risk of misdiagnosis when using diagnostic instruments CT or MRI, where the presence of steel implant in that area will result in poor display.

With each surgery intended to remedy consequences of multiple injuries, severe complications may occur. These complications include, inter alia, genitourinary, gastrointestinal and vascular disorders including thrombosis; bronchopulmonary disorders including embolism, myocardial infarction and death.

Reprocessing limitation

Implants are disposables for humane use. If the implants were prepared for surgery and not contaminated, they can be re-prepared for use by cleaning, disinfection and sterilization process as described below.

INSTRUCTIONS

Place of use:

Professional workplace of health care providers, namely surgical, trauma and orthopedic departments.

Storage and transport:

Prevent deformation and scratches to the implant; avoid contact with other metallic materials and chemicals when handling the implant.

Preparation for cleaning

Implants are cleaned by the manufacturer, and therefore need not be cleaned or disinfected. Should it be necessary to re-sanitize an unused implant, use a disinfectant that contains no chloride ions and is suitable for steel instruments. The implant must not be damaged (scratched) during disinfection.

Cleaning – automated

No special requirements.

Cleaning – manual

No special requirements.

Disinfection:

We recommend either the equipment for washing and thermal or thermochemical disinfection or hand washing of instruments after chemical disinfection with virucidal agent, while the two stages can be combined using disinfectants with combined washing effect. Disinfectant can be used in accordance with the instructions on the label of the product. In the preparation of disinfectants and cleaning solutions, the procedures specified by the manufacturer must be followed. Implants that are designed for disassembly are disassembled and each part is further considered as separate. After disinfection with implants rinse with demineralized water to remove chemical residues remaining on the implant.

Recommended products: Sekusept pulver, NeodisherseptoMED

Drying:

If drying is part of cycle for washing and disinfecting device, a temperature of 120 °C should not be exceeded.

Maintenance:

It is not performed.

Inspection and testing of function:

If cleaning and disinfection are applied before use, the implant should be inspected visually for damage.

Packing:

It is possible to use standard packaging. Ensure that the container is large enough to accommodate the implant without stressing the seals. The implants are placed in containers or individually in sterilization trays, which are then packaged (cutting edges and tips must be protected).

Sterilization:

We recommend that implants are sterilized separately from other instruments. Unless otherwise stated in a special notification, it is advisable to use steam sterilization in a steam sterilizer equipped with antibacterial filter at 121 °C, pressure 205 kPa for 20 minutes or at 134 °C, pressure 304 kPa for 10 minutes. If the special notification includes recommendations on hot air sterilization, it must be carried out in devices with forced air circulation at 160 °C for 60 minutes or at 170 °C for 30 minutes or at 180 °C for 20 minutes, unless otherwise indicated. Hot air sterilizer should be opened after the end of the cycle after cooling to 80 °C at least.

Storage:

The implant must be stored in a dry place without direct exposure to chemicals.

Other information:

When sterilizing multiple implants in one autoclave cycle, be sure not to exceed the maximum load of the sterilizer. Cartridges, boxes, screens and stands can be used to store, transfer and sterilize medical devices. Unpunched cartridges, boxes and screens cannot be used as a packaging material (container) used to sterilize medical devices. Open those before sterilization.

Implant disposal:

The implant is considered hazardous waste after use. User is responsible for implementing measures of safe handling and disposal of the product. After drying, used implants are discarded and disposed of as potentially hazardous waste.

The joint-stock company MEDIN validated the above mentioned methods of cleaning, disinfection, drying and sterilization.Through the validation were all procedures recognized as QUALIFIED for the preparation on medical devices for repeated using. The user is responsible for securing of necessary materials, equipment and staff training in accordance with these specified procedures to reach the required effect. The user can use some other methods but he must consider and judge (evaluate) the effectivity of these procedures considering the potential unwanted impacts.

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