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<div data-bbox="347 6

- In the patient do not use implants made of various materials because a contact between various types of metals can accelerate the corrosion process which may lead to premature decrease in efficiency and lifetime of the implant; it may lead to the increased amount of metal compounds that released in the patient's body and may cause a toxic reaction.
- At the medical facility keep available all shapes and sizes of the implants that are likely to be used.
- Take into account additional preventive measures that reduce the weight load to minimize the premature loading of the means used for fracture fixation until the bone connection is not visible on the X-ray.
- Decide on the removal of the implant for internal fixation, taking into account the risks for the patient due to the repeated surgery, difficulties during and after the removal of the implant and the benefits of the implant removal for the patient.
- The patients who smoke or who usually use product containing nicotine have a higher number of unhealed fractures. These patients should be warned against the possible consequences. The incidence of implant failure is more likely to occur among heavy, physically active people.
- The manufacturer does not bear liability if the product is used differently than described in the manual, eventually surgery procedure.
- The manufacturer does not bear liability if the product is used in combination with implant so of other manufacturers.

#### Instruction for patients:

- Tell the patient to precisely follow the instructions of his/her attending physician.
- Inform the patient about the surgery, possible risks and rehabilitation programmes.
- Inform the patient when he/she can fully load the bone with an internal fixation to prevent the damage of the implant.
- Tell the patient that the bone during the treatment is susceptible to damage, especially in places where the implant ends.
- Inform the patient about the possible removal of the implant.
- Inform the patient to retain a reasonable level of activities after the removal of the implant, because the bone may fracture again until the bone is not fully healed and until the bone cavity is filled.

#### Possible side effects:

- Not finished connection of the fracture ends, osteoporosis, suppressed revascularization and poor bone formation can cause a loosening, bending, breakage or premature loss of the strong bone fixation.
- Poor connection of the fracture resulting from improper alignment
- Increased response of the fibrous tissue surrounding the fracture due to unstable fragmented fractures
- Early or late infection, both depth and surface infection.
- Avascular necrosis.
- Shortening of the affected side of the bone (fracture), in children on the contrary the outgrow of the bone length and thus the limb
- Subclinical damage of nerves may occur as a result of a surgical trauma.
- Allergies to nickel
- Risk of misdiagnosis when using the diagnostic instruments CT or MRI, where the presence of the steel implant in the area may cause poor display.

Severe complications may occur during or after the surgery to eliminate the consequences of multiple injuries. These complications include, inter alia, genitourinary, gastrointestinal and vascular disorders, including thrombosis, bronchopulmonary disorders, including embolism, myocardial infarction and death.

#### Restriction of reoccurring processing

The implants are intended for a single use in humans. If the implants were prepared for surgery and are not contaminated, they can be ready for use again using the process of cleaning, disinfection and sterilization described below

#### INSTRUCTIONS

##### WARNING:

The implant is supplied NON-STERILE. It is packaged in a polyethylene bag. It comes cleaned and disinfected from the manufacturer. Before first use, remove the implant from the package and perform sterilization as prescribed. If the package is broken, it is necessary to return the implant to the manufacturer for repackaging.

##### Place of use:

Professional workplace of health care providers, namely surgical, trauma and orthopedic departments.

##### Storage and transport:

Prevent deformation and scratches to the implant; avoid contact with other metallic materials and chemicals when handling the implant.

##### Preparation for cleaning

Implants are cleaned by the manufacturer, and therefore need not be cleaned or disinfected. Should it be necessary to re-sanitize an unused implant, use a disinfectant that contains no chloride ions and is suitable for steel instruments. The implant must not be damaged (scratched) during disinfection.

##### Disinfection:

We recommend either the equipment for washing and thermal or thermochemical disinfection or hand washing of instruments after chemical disinfection with virucidal agent, while the two stages can be combined using disinfectants with combined washing effect. Disinfectant can be used in accordance with the instructions on the label of the product. In the preparation of disinfectants and cleaning solutions, the procedures specified by the manufacturer must be followed. Implants that are designed for disassembly are disassembled and each part is further considered as separate. After disinfection with implants rinse with demineralized water to remove chemical residues remaining on the implant.

*Recommended products: Sekusept pulver, NeodisherseptoMED*

##### Drying:

If drying is part of cycle for washing and disinfecting device, a temperature of 120 °C should not be exceeded.

##### Maintenance:

It is not performed.

##### Inspection and testing of function:

Before use, implant's condition should be inspected visually.

If cleaning and disinfection are applied before use, the implant should be inspected visually for damage.

##### Sterilization:

We recommend that implants are sterilized separately from other instruments.

Unless otherwise stated in a special notification, it is advisable to use steam sterilization in a steam sterilizer equipped with antibacterial filter at 121 °C, pressure 205 kPa for 20 minutes or at 134 °C, pressure 304 kPa for 10 minutes.

When sterilizing multiple implants in one autoclave cycle, be sure not to exceed the maximum load of the sterilizer. Cartridges, boxes, screens and stands can be used to store, transfer and sterilize medical devices. Unpunched cartridges, boxes and screens cannot be used as a packaging material (container) for sterilization of medical devices. Open those before sterilization.

##### Storage:

The implant must be stored in a dry place without direct exposure to chemicals.

##### Other information:

The joint-stock company MEDIN validated the above mentioned methods of cleaning, disinfection, drying and sterilization. Through the validation were all procedures recognized as QUALIFIED for the preparation on medical devices for repeated using. The user is responsible for securing of necessary materials, equipment and staff training in accordance with these specified procedures to reach the required effect. The user can use some other methods but he must consider and judge (evaluate) the effectiveness of these procedures considering the potential unwanted impacts.

##### Implant disposal:

Used implant is considered hazardous waste. User is responsible for implementing measures of safe handling and disposal of the product.

**After drying, used implants are discarded and disposed of as potentially hazardous waste.**

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