VAROVÁNÍ

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
END CAPS

Description of the implant

End caps are made of titanium alloy Ti6Al4V (ISO 5832-3).


The end cap has a cylindrical shape with M8 external thread. There is inner hexagon 5 mm on the end cap, through which the end cap introduced into the nail with a screwdriver.

End cap for femoral retrograde nail

NAME	STEEL	TITANIUM
End cap M8xL5 mm	–	129 78 4143

End cap for tibiototalcaneal arthrodesis retrograde nail

NAME	STEEL	TITANIUM
End cap M8xL10 mm	–	129 70 6753

Compatibility

End caps are designed for use with the following nails only. Special instrumentation is used for their insertion

Nails:

Retrograde femoral nail (PLO222)

397 129 78 6793 – 397 129 78 6883

Lengths **175–350 mm** (175/200/250/300/350)

Tibiototalcaneal arthrodesis retrograde nail (PLO223)

397 129 70 0013 – 397 129 70 0093

Lengths **180–250 mm** (180/220/250)

397 129 70 2903 – 397 129 70 2983

Lengths **180–250 mm** (180/220/250)

Instrumentation:

Retrograde femoral nail (PLO234)

397 139 09 0575 – Instrumentation for retrograde femoral nails with sieve

Tibiototalcaneal arthrodesis retrograde nail (PLO231)

397 139 09 0595 – Instrumentation for tibiototalcaneal arthrodesis retrograde nail with sieve

Indication

The end cap is designed to close the end of the nail and lock the last screw at the same time.

Contraindication

Are not.

Instructions for use

Insert the end cap into the distal end of the nail with a screwdriver. Tightening will lock the last screw in the nail.

A detailed description of the implantation is provided in the **Operation procedure – Retrograde femoral nail, Retrograde nail for tibiototalcaneal arthrodesis.**

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

WARNING

Implants (nails, screws, end caps) from different materials may not be combined

STEEL x TITANIUM ❌

TITANIUM x TITANIUM or STEEL x STEEL ✅

Products may not be combined with products by other manufacturers.

Reprocessing limitation

Implants are disposables for humane use. If the implants were prepared for surgery and not contaminated and mechanically damaged, they can be re-prepared for use by cleaning, disinfection and sterilization process as described below.

INSTRUCTIONS

WARNING: The implant is supplied NON-STERILE. It is packaged in a polyethylene bag. It comes cleaned and disinfected from the manufacturer. Before first use, remove the implant from the package and perform sterilization as prescribed. If the package is broken, it is necessary to return the implant to the manufacturer for repackaging.

Place of use: Professional workplace of health care providers, namely surgical, trauma and orthopedic departments.

Storage and transport: Prevent deformation and scratches to the implant; avoid contact with other metallic materials and chemicals when handling the implant.

Preparation for cleaning: Implants are cleaned by the manufacturer, and therefore need not be cleaned or disinfected. Should it be necessary to re-sanitize an unused implant, use a disinfectant that contains no chloride ions and is suitable for steel instruments. The implant must not be damaged (scratched) during disinfection.

Disinfection: We recommend either the equipment for washing and thermal or thermochemical disinfection or hand washing of instruments after chemical disinfection with virucidal agent, while the two stages can be combined using disinfectants with combined washing effect. Disinfectant can be used in accordance with the instructions on the label of the product. In the preparation of disinfectants and cleaning solutions, the procedures specified by the manufacturer must be followed. Implants that are designed for disassembly are disassembled and each part is further considered as separate. After disinfection with implants rinse with demineralized water to remove chemical residues remaining on the implant. *Recommended products:* Sekusept pulver, NeodisherseptoMED

Drying: If drying is part of cycle for washing and disinfecting device, a temperature of 120 °C should not be exceeded.

Maintenance: It is not performed.

Inspection and testing of function: Before use, implant's condition should be inspected visually. If cleaning and disinfection are applied before use, the implant should be inspected visually for damage.

Sterilization: We recommend that implants are sterilized separately from other instruments. Unless otherwise stated in a special notification, it is advisable to use steam sterilization in a steam sterilizer equipped with antibacterial filter at 121 °C, pressure 205 kPa for 20 minutes or at 134 °C, pressure 304 kPa for 10 minutes. When sterilizing multiple implants in one autoclave cycle, be sure not to exceed the maximum load of the sterilizer. Cartridges, boxes, screens and stands can be used to store, transfer and sterilize medical devices. Unpunched cartridges, boxes and screens cannot be used as a packaging material (container) for sterilization of medical devices. Open those before sterilization.

Storage: The implant must be stored in a dry place without direct exposure to chemicals.

Other information: The joint-stock company MEDIN validated the above mentioned methods of cleaning, disinfection, drying and sterilization. Through the validation were all procedures recognized as QUALIFIED for the preparation on medical devices for repeated using. The user is responsible for securing of necessary materials, equipment and staff training in accordance with these specified procedures to reach the required effect. The user can use some other methods but he must consider and judge (evaluate) the effectivity of these procedures considering the potential unwanted impacts.

Implant disposal: Used implant is considered hazardous waste. User is responsible for implementing measures of safe handling and disposal of the product.

After drying, used implants are discarded and disposed of as potentially hazardous waste.

Manufacturer: CZ43378030 Number and date of last revision: R01/2014-07-02
MEDIN, a.s. | Vlachovická 619 | CZ 592 31 Nové Město na Moravě | Czech Republic
tel.: +420 566 684 336 | fax: +420 566 684 385 | e-mail: export@medin.cz

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